IntraTemp™

Solution Warmer

OPERATOR'S MANUAL

IntraTemp™ Solution Warmer Operator's Manual

TABLE OF CONTENTS

1.	Introduction	3
	• Intended Use / Indications for Use	
	• Contraindications	
	• Important Labels and Indicators	
	• WARNINGS AND PRECAUTIONS	
2.	Specifications	4
3.	Operating Characteristics	7
	• Power Connection	
	• Heating Plate	
	• Fluid Basin-ThermaBasin™ Technology	
	• Heater Safety Switch	
	• Temperature Control System	
	• Upper Basin Indicator Light and Digital Display	
	• Heated Bottle Pocket™ Technology	
	• Drain Tank	
4.	Unpacking Instructions	8
	• Unpacking the IntraTemp Solution Warmer	
5.	Installation Instructions	8
	• Tools Required	
	• Connecting the Power Cord	
	• Basin Heater Controller Setup	
	• Drain Tank Installation	
6.	Operating Instructions	11
	• Steps to use the IntraTemp Solution Warmer	
	• Installing a ThermaBasin	
7.	Troubleshooting	12
8.	Maintenance	12
	• Care of IntraTemp Solution Warmer After Each Use	
	• Bi-Weekly Maintenance of IntraTemp Solution Warmer	
	• Fuse Replacement Instructions	
	• IntraTemp Basin Temperature Verification	
9.	Contact Information	13

Rx ONLY CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

1. INTRODUCTION

- 1.1 **INTENDED USE / INDICATIONS FOR USE** The IntraTemp™ Solution Warmer (ITW6L) is designed to warm and maintain the temperature of surgical solutions prior to their use.
- 1.2 **CONTRAINDICATIONS:** The IntraTemp Solution Warmer is NOT intended for use as a blood warmer, as a tissue bath, or as a patient (e.g., an infant) bath.
- 1.3 **IMPORTANT LABELS AND INDICATORS**
 - Digital Temperature Display and Green Power-On Light beside Basin.
 - Digital Temperature Displays above each Bottle Pocket.
 - Power Switch- “I” ON, “0” OFF.
 - Serial Number & Ratings found beside power switch on lower rear of the IntraTemp cabinet.
 - **Explanation of Labels:**

Symbol	Meaning	Application
	Electric Shock Hazard	The front panel should not be removed. No user serviceable parts inside.
	Hot Surface	Burn Hazard. Allow this surface to cool before cleaning.
No Direct Patient Contact with Heating Unit!	No Patient Should Touch The Unit!	Direct patient contact is contraindicated! The IntraTemp™ is not intended as a blood warmer, tissue bath, or patient bath.
	Protective Earth	Source of Earth Ground for the device.
	Caution	Caution Sign Adjacent to Basin Heater and Bottle Pocket- further warning that the heater surface may be hot. Caution Sign Adjacent to Rating Label - use only specified fuse.
	3rd Party IEC 60601-1 3rd Edition Certification Listing	Located on Rating Label - Certified that IntraTemp conforms to IEC 60601-1 3rd Edition Standard as tested by a 3rd Party Laboratory.
	CE Mark	Located on Rating Label - Indicates that IntraTemp conforms to European Union Directives
	Symbol for marking of Electrical and Electronic Equipment (EEE).	Located on Rating Label - Disposal of Old Electrical & Electronic Equipment (Applicable in the European Union and other European countries with separate collection systems) This symbol on the product or on its packaging indicates that this product shall not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product. ECOLAB has contracted with one or more recycling/disposal companies in each EU Member State, Please contact your ECOLAB contact for recycling/disposal information.
	Manufacturer Date of Manufacture	Located on Rating Label

Symbol	Meaning	Application
	Authorized Representative in the European Community	Located on Rating Label
	Catalog Number	Located on Rating Label
	Serial Number	Located on Rating Label
	Follow Instructions for Use	Located on Rating Label
	Fuse	Located on Rating Label
	Loaded Mass	Located on Rating Label – Indicates equipment mass including working load
	Do not push with casters locked or blocked to prevent tipping	Located on Rating Label

1.4 WARNINGS AND PRECAUTIONS

- **NO DIRECT PATIENT CONTACT! DO NOT PLACE A PATIENT** (e.g., an infant) OR ANY PART OF A PATIENT INTO, ON, OR AGAINST THE IntraTemp OR INTO DIRECT CONTACT WITH THE WARMING PLATE LOCATED IN THE UPPER BASIN RECEPTACLE. PROLONGED CONTACT WITH HEATING SYSTEMS MAY CAUSE THERMAL INJURY.
- IF THE WARMER IS POWERED ON WITH A THERMABASIN INSTALLED, THE BASIN HEATING PLATE WILL HEAT UP EVEN IF NO FLUID IS POURED INTO THE BASIN. WITHOUT FLUID IN THE BASIN, THE THERMABASIN HEATING PLATE MAY BECOME HOT ENOUGH TO CAUSE THERMAL INJURY. DO NOT OPERATE WITHOUT A MINIMUM OF 1.5 LITERS OF SOLUTION IN THE THERMABASIN.
- THE BOTTLE POCKET HEATED ACCESS SLEEVES ARE ALWAYS HEATING WHEN POWER IS ON. CARE SHOULD BE TAKEN TO NOT TOUCH THE INTERIOR METAL SURFACES AS THEY CAN BE HOT. OBJECTS OTHER THAN SALINE BOTTLES OR BAGS SHOULD NEVER BE PLACED IN THE POCKETS
- POSITION INTRATEMP TO ALLOW EASY REMOVAL OF POWER CORD FROM POWER OUTLET. POWER CORD MUST BE DISCONNECTED FROM OUTLET TO ISOLATE INTRATEMP FROM EXTERNAL POWER SUPPLY.
- ELECTRIC SHOCK HAZARD- DO **NOT** REMOVE PROTECTIVE PANELS. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL.
- RISK OF FIRE. REPLACE FUSE AS MARKED:
- CANADA AND US: GROUNDING RELIABILITY CAN ONLY BE ACHIEVED WHEN THIS EQUIPMENT IS CONNECTED TO AN EQUIVALENT RECEPTACLE MARKED “HOSPITAL ONLY” OR “HOSPITAL GRADE”
- INTRATEMP IS DESIGNED FOR USE BY TRAINED HEALTHCARE PROFESSIONALS.
- INTRATEMP IS DESIGNED FOR USE IN INDOOR AMBIENT TEMPERATURE AND HUMIDITY CONDITIONS. (10-40 °C; 30-75 %RH).
- WARNING: DO NOT TRANSPORT FULLY LOADED AT AN INCLINE GREATER THAN 5°.
- INTRATEMP IS NOT DESIGNED FOR USE IN AN OXYGEN RICH ENVIRONMENT.
- INTRATEMP IS NOT DESIGNED FOR USE IN THE PRESENCE OF FLAMMABLE ANESTHETICS.
- WARNING: NO MODIFICATION OF THIS EQUIPMENT IS ALLOWED.

2. SPECIFICATIONS

- Power: 120 Volt AC, 60Hz, 5 Amps/230V AC, 50Hz, 2.6 Amps.
- Voltage Sense Circuit Software Revision: ITW6L-_-MSIBOARD.X.PRODUCTION.hex

- North American fuse: F 6.3AL, 250 VAC (5 x 20 mm).
- European fuse: F 3.15AL, 250 VAC (5 x 20 mm).
- Temperature Scale: Celsius or Fahrenheit selectable by qualified service personnel.
- Basin Temperature Set Point: selectable by qualified service personnel in a range from 37 to 43 ± 1°C (98 to 110 ± 1.5°F).
- Heated Bottle Pocket Receptacles Temperature Set-Point: Factory setting per Customer Specification up to 43 ± 1°C (110 ± 2°F).
- Use only with Ecolab ThermaBasin; product code ITB100.
- Minimum fluid volume: 1.5 Liters.
- Equipment mass including working load is 43 Kg
- EMC Compliance: The ITW6L has been evaluated with respect to electric shock, fire and mechanical hazards in accordance with EN IEC 60601-1, AAMI ES60601-1, CAN/CSA C22.2 – No. 60601-1.

Guidance and manufacturer's declaration – electromagnetic emissions		
The ITW6L is intended for use in the electromagnetic environment specified below. The customer or the user of the ITW6L should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The ITW6L uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The ITW6L is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 6100-3-2	Class B	
Voltage fluctuations flicker emissions IEC 6100-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The ITW6L is intended for use in the electromagnetic environment specified below. The customer or the user of the ITW6L should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 6100-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 6100-4-4	±2 kV for power supply lines ±2 kV for input/output lines	±2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 4100-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.

Guidance and manufacturer's declaration – electromagnetic immunity

The ITW6L is intended for use in the electromagnetic environment specified below. The customer or the user of the ITW6L should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC6100-4-11	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ITW6L requires continued operation during power mains interruptions, it is recommended that the ITW6L be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 6100-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. Portable and mobile RF communications equipment should be used no closer to any part of the ITW6L, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF IEC 6100-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d=1.2\sqrt{P}$
Radiated RF IEC 6100-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE— U_T is the a.c. mains voltage prior to application of the test level.

NOTE 1—At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ITW6L is used exceeds the applicable RF compliance level above, the ITW6L should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ITW6L.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between
portable and mobile RF communications equipment and the ITW6L**

The ITW6L is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ITW6L can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ITW6L as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1—At 80 MHz to 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Catalog Number	Description
ITW6L	IntraTemp Solution Warmer
ITNA15	Power Cord Kit, IntraTemp, North American Hospital Grade, 4.6 m (15 ft)
ITNA25	Power Cord Kit, IntraTemp, North American Hospital Grade, 7.6 m (25 ft)
ITUK15	Power Cord Kit, IntraTemp, United Kingdom, 4.6 m
ITUK25	Power Cord Kit, IntraTemp, United Kingdom, 7.6 m
ITEU15	Power Cord Kit, IntraTemp, European Union, 4.6 m
ITEU25	Power Cord Kit, IntraTemp, European Union, 7.6 m
ITB100	Ecolab ThermaBasin: Basin with Attached Drape for use with the IntraTemp Solution Warmer

3. OPERATING CHARACTERISTICS

3.1 Power Connection

The IntraTemp Solution Warmer is designed to work when connected to a 120 VAC 60 Hz Power Source (North America) or a 230 VAC 50 Hz Power Source (UK, Continental Europe). A power cord specific to an applicable region must be ordered with the IntraTemp Solution Warmer.

3.2 Heating Plate

The IntraTemp Solution Warmer warms fluid contained in a basin via a specially designed ThermaBasin heating plate. The heating plate is located at the top of the warmer in the basin receptacle. After use, allow the heating plate to cool before cleaning. Avoid contact until the plate cools.

3.3 Fluid Basin - ThermaBasin™ Technology

IntraTemp is designed to operate with a proprietary Ecolab ThermaBasin with a unique thermal well molded into its bottom. This molded thermal well protrudes upward from the disposable ThermaBasin bottom and mates with the warmer's Safety Switch thermocouple to measure and maintain fluid temperature. Only Ecolab ThermaBasin can be used with the IntraTemp warmer. The use of any other basin may damage or disable the IntraTemp warmer.

3.4 Heater Safety Switch

The ThermaBasin heating plate Safety Switch Thermocouple has three positions:

- Fully up- no ThermaBasin in place - no heat
- Middle position- ThermaBasin in place - heat conducted to ThermaBasin
- Fully down - no heat

Care should be exercised that the ThermaBasin is properly located in the basin receptacle with the basin's thermal well placed securely over the warmer's Safety Switch thermocouple.

3.5 Temperature Control System

To ensure the safe use of warm fluids the basin heater controls may be set at a temperature between 37° and 43°C (98° - 110°F) and the temperature scale may be set to either Celsius or Fahrenheit by qualified service personnel. Access to the control buttons by the user is restricted by a controller cover.

3.6 Upper Basin Indicator Light and Digital Display

The green "Power On" indicator light glows to indicate that power is available to the upper basin heater and that a ThermaBasin is properly installed for use.

The digital display on the top surface of the IntraTemp displays the temperature of the solution in the fluid basin.

If the ThermaBasin is not properly seated, neither the green indicator light nor the digital display will be activated; the basin heating system is not actuated.

3.7 Heated Bottle Pocket Technology

The side housing of the solution warmer has two (2) Bottle Pocket heated access sleeves to hold standard 1 or 1.5 liter sterile saline solution bottles or 1 liter bags.

The Bottle Pockets cannot control a fluid bottle or bag to an exact temperature; they are designed so that when a fluid container is in place it will be maintained at or very near a steady temperature.

3.8 Drain Tank

A Drain Tank is provided to capture any fluid spilled into the upper housing receptacle. A tube that runs inside the unit connects the drain hole located in the upper housing receptacle to the drain tank.

4. UNPACKING INSTRUCTIONS

4.1 Unpacking the IntraTemp Solution Warmer:

1. DO NOT CUT THE BOX. Cut the strapping bands and remove.
2. Cut the tape holding the top of the box closed. Lift the box flaps.
3. Remove the Drain Tank from the top of the packaging material and set aside for installation onto the unit.
4. Remove the packaging material from the top of the unit.
5. Lift the box off the IntraTemp machine.
6. Remove IntraTemp from the lower packaging.
7. Save the shipping material if a unit is to be returned.

NOTE: Locate the Power Cord Kit supplied in a separate box.

5. INSTALLATION INSTRUCTIONS

NOTE: Procedures described in this section should only be performed by qualified service personnel or Ecolab/Microtek Factory Representative

5.1 Tools Required:

The procedures described in this section will require the following tools not included in the IntraTemp Solution Warmer packaging:

- #1 Phillips Head Screw Driver
- #2 Phillips Head Screw Driver

- Flat Head Screw Driver

5.2 Installing Fuse:

1. Remove the fuses from the power cord kit.
2. Install the appropriate fuse for the voltage source used.
 - North American fuse: F 6.3AL, 250 VAC (5 x 20 mm).
 - European fuse: F 3.15AL, 250 VAC (5 x 20 mm).
3. Use a small flathead screw driver to depress the latch and open the fuse drawer as shown in Figure 1.

Figure 1: Fuse drawer opening

4. Install a fuse into the rear clip as shown in Figure 2.
5. Install a second fuse into the spare fuse holder located at the front of the drawer as shown in Figure 2.
6. Close fuse drawer and ensure the drawer is securely latched.

Figure 2: Fuses installed

5.3 Connecting the Power Cord:

1. Remove the power cord from the power cord kit.
2. Check to ensure that the plug on the power cord is compatible with the applicable power receptacle in the region of use. If the power cord is not compatible contact an Ecolab/Microtek representative.
3. Attach the power cord to the power entry module located below the main power switch on the back of the IntraTemp Solution Warmer.
4. Remove the power cord retention bracket and two mounting screws from the power cord kit.
5. Using a #1 Phillips Head Screw Driver and furnished screws connect the power cord retention bracket to the IntraTemp Solution Warmer via the threaded connectors located to each side of the power inlet module (see Figure 3).

Figure 3: Power Cord Installation

5.4 Basin Heater Controller Setup:

To access the basin warmer controller use a #2 Phillips Head Screw Driver to remove the two screws securing the controller cover to the IntraTemp main body (see Figure 4).

Figure 4: Controller Cover Removal

Figure 5: Temperature Adjustment

Plug the power cord into a power receptacle.

Turn on power switch.

The controller indicates the basin warmer current temperature of the basin in the upper display and the set-point temperature in the lower display.

To change the temperature scale on all of the temperature displays perform the following step:

- Press the EZ1 button on the Watlow Controller to toggle between Celsius and Fahrenheit.

To adjust the basin warmer controller set-point between 37 and 43°C (98° - 110°F):

- Use the appropriate UP ▲ or Down ▼ buttons to adjust the basin warmer controller set point as shown in Figure 3.

Visually check the Bottle Pocket Warmer Display to ensure that a temperature near room ambient is displayed in the preferred temperature scale.

Turn off the power switch and remove the power plug from the wall receptacle.

Reattach the controller cover using a #2 Phillips Head Screw Driver and the two screws supplied with the unit.

5.5 Drain Tank Installation:

Remove packaging from the supplied Drain Tank.

Grasp the Drain Tank utilizing the furnished recess on the bottom of the tank, and place on the rails located on the bottom of the IntraTemp Solution Warmer.

NOTE: Insert the drain tank with the intake notch side first as shown in Figure 6.

Insert the drain tank until it is fully seated; positioned under the warmer without overhang as shown in Figure 7.

Figure 6: Drain Tank Installation

Figure 7: Drain Tank fully seated

6. OPERATING INSTRUCTIONS

6.1 Follow these simple steps to use the IntraTemp Solution Warmer:

1. Plug the power cord into a power receptacle. It is recommended to ensure sufficient capacity is available on the circuit being used.
2. Turn on power switch. This actuates the Bottle Pocket heaters automatically; this also allows the user to actuate the basin heater when ThermaBasin is correctly installed.
3. Correctly install a ThermaBasin (See section 6.2 below).
4. Pour a minimum of 1.5 L of solution into the ThermaBasin.
5. The green Power-On indicator light, digital fluid temperature display and basin heater will activate when a ThermaBasin is properly installed.
6. Saline bottles or bags may be placed into the Bottle Pockets as needed.
7. At the end of a procedure, removing the ThermaBasin will turn off power to the basin heater, shutting off the green power indicator light and basin temperature display.
8. Turn off the unit's power switch.
9. Discard ThermaBasin per hospital guidelines.
10. After use, allow the heating plate to cool before cleaning.
11. Clean and dry the ThermaBasin receptacle after each use. Do not use abrasive cleaners or cloths. It is particularly important to avoid scratching the ThermaBasin receptacle and that it be clean and dry between uses.

NOTE: A small amount of water can accumulate in bottom of the basin receptacle due to condensation during use or due to fluid spilled when removing the ThermaBasin after use.

6.2 Installing a ThermaBasin

1. Use only Ecolab ThermaBasin with the integrated thermal well and drape. The IntraTemp Solution Warmer functions only with Ecolab ThermaBasin.
2. Remove a ThermaBasin with attached drape from the sterile package per accepted aseptic technique.
3. With the drape collected on top of the ThermaBasin and the bottom of the ThermaBasin exposed, place the ThermaBasin into the basin receptacle on top of the warmer.
4. Align the basin's thermal well with the Safety Switch in the bottom of the warmer's basin receptacle. Match the basin's locator slot with the locator ridge in the basin receptacle. Press the basin bottom and the basin's thermal well onto the solution warmer's Safety Switch. Pour fluid into basin.
5. When the basin is properly installed both the green Power On indicator light and the digital fluid temperature display will illuminate beside the basin this verifies that the basin heater is actuated.
6. The heater will not actuate until the basin is properly seated and the warmer's power switch is in the ON position. Firmly depress the basin's thermal well onto the Safety Switch. If the basin is properly installed and the warmer's power switch is on but the power indicator light and fluid temperature-display do not illuminate, please unplug the unit and contact Ecolab/Microtek Customer Service.

Note: If added fluid is significantly above or below the temperature of steady state fluid already in use in the ThermaBasin, it may take a short period of time for the system to correctly measure and display the fluid temperature.

7. TROUBLESHOOTING

If these steps do not correct the problem, please contact Ecolab/Microtek Customer Service.

PROBLEM	Indicator light and digital fluid temperature display do not illuminate when power is on.
SOLUTION	Check ThermaBasin is properly seated. Check power cord. Check fuse.
PROBLEM	Fluid is not being heated.
SOLUTION	Ensure that an Ecolab ThermaBasin is being used. Check that basin is properly seated over the Safety Switch. Both the green Power On indicator light and the digital fluid temperature display beside the basin should be illuminated.
PROBLEM	Fluid does not reach set temperature or has unusual fluctuation reading.
SOLUTION	Check that basin is properly seated and the basin bottom is in direct contact with heating plate. Make sure that drape material is not between bottom of basin and heating plate.

8. MAINTENANCE

8.1 Care of IntraTemp Solution Warmer after Each Use:

The IntraTemp Solution Warmer is designed to require minimal maintenance. Clean and dry the ThermaBasin receptacle after each use. Do not use abrasive cleaners or cloths. It is particularly important to avoid scratching the ThermaBasin receptacle and that it be clean and dry between uses.

8.2 Bi-Weekly Maintenance of IntraTemp Solution Warmer:

The IntraTemp basin receptacle should be scrubbed with a non-abrasive cleanser every two weeks. This prevents the build-up of corrosion or staining in the basin receptacle which may be caused by normal condensation during use.

The power cord should be inspected for any unusual sign of wear or abuse. If abuse or wear are noted, please contact Ecolab/Microtek Customer Service to have the power cord replaced.

8.3 Fuse Replacement Instructions:

Unplug the power cord from the power receptacle. The following steps are required to replace the main power fuse:

- Use a #1 Phillips Head Screw Driver to remove the two screws securing the power cord retention bracket to the IntraTemp ITW6L, and set aside.
- Unplug the power cord from the power inlet module.
- Use a small flathead screw driver to depress the latch on the fuse drawer as shown in Figure 8.
- Remove the rear fuse from the clip shown in Figure 9.
- Replace the fuse with the spare fuse located at the front of the fuse drawer. REPLACE FUSE AS MARKED.
- Close fuse drawer and ensure the drawer is securely latched.
- Plug the power cord into the power inlet module.

Figure 8: Fuse Drawer Opening

Figure 9: Fuse Drawer Opening

- Using a #1 Phillips Head Screw Driver and furnished screws connect the power cord retention bracket to the IntraTemp Solution Warmer via the threaded connectors located to each side of the power inlet module.

8.4 IntraTemp Basin Temperature Verification:

The IntraTemp basin solution temperature should be verified every five (5) years. The following items are required to perform this verification:

- Calibrated Fluke 51II Digital Thermometer or equivalent.
- Shrouded 6" long Type T Thermocouple Probe.

Install a ThermaBasin. Add two (2) liters of saline or water into the ThermaBasin. Turn the unit on and allow the fluid to heat for two (2) hours.

Place the tip of the Thermocouple Probe into the fluid (without contacting the basin) at a point ~1" (2.5 cm) above the center of the basin. The temperature indicated by the calibrated thermometer should be at Set-Point $\pm 1^{\circ}\text{C}$ ($\pm 1.5^{\circ}\text{F}$). If the temperature is not within the specified range, contact Ecolab/Microtek Customer Service.

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© 2014 ECOLAB INC. o sus afiliados. Todos los derechos reservados
Medical, Inc. 602
Columbus, MS 39702 U.S.A.
Microtek Medical B.V.