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Operating instructions

Dermados Station Touchfree

Automatic “touchfree” dosing dispenser for hand sanitizing agent

Dermados Station
MAN048961, Rev. 2-09.2020
16.09.2020

ENGLISH

Table of contents

1	General information	3
1.1	References to Instructions of Use.....	3
1.2	Transportation.....	6
1.3	Repairs / returns to Ecolab Engineering GmbH.....	7
1.4	Packaging.....	8
1.5	Storage.....	9
1.6	Equipment marking – Rating plate.....	9
1.7	Warranty.....	10
1.8	Contacts.....	10
2	Safety	12
2.1	General safety advice.....	12
2.2	Use for the intended purpose.....	12
2.3	Arbitrary conversion and spare parts production.....	16
3	Scope of delivery	17
4	Functional description	18
4.1	Variants.....	19
5	Installation and assembly	21
5.1	Setup.....	22
5.2	Assembly.....	23
6	Commissioning	25
7	Operation - Hand sanitizing	28
8	Description and rectification of malfunctions	30
9	Service	31
9.1	Service.....	32
10	Spare parts and accessories	38
10.1	Dosing container.....	39
11	Conversion - canister to bottle or reactivation of canister version	40
12	Technical specifications	42
12.1	Type plate.....	42
12.2	Dimensions.....	43
13	Decommissioning, disassembly, environmental protection	45
14	EC Declaration/Declaration of Conformity	46
	Appendix	47
A	Safety data sheets.....	48

1 General information

1.1 References to Instructions of Use

CAUTION!

Observe the Instructions of Use!

Please read through these Instructions of Use prior to putting the dispenser into use and observe the device safety information and the information on the intended use. Keep the Instructions of Use so that you can quickly find and rectify any malfunctions whose cause you can find using the information (↪ *Chapter 8 'Description and rectification of malfunctions' on page 30*). We recommend these Instructions of Use be reviewed again from time to time so as not to forget routine jobs and to keep that expertise refreshed.

The German version of the are the **Instructions of Use in the original language**, that are legally relevant. **All the other languages are translations.**

- Illustrations in this manual are provided for basic understanding and may deviate from the actual product.
- Prior to assembly, putting into first use and doing any maintenance or repair work, it is important that the relevant chapters of these Instructions of Use be read, understand and observed.

Complete Instructions of Use for downloading

The latest and complete Instructions of Use are available online.

To download the instructions with a PC, tablet or Smartphone, use the link below or scan the QR code shown.

This QR code is also located on the back of the Dermados Station!

Downloading the Instructions of Use:

https://www.ecolab-engineering.de/fileadmin/download/bedienungsanleitungen/institutional/MAN048961_Dermados_Station.pdf

Videolink

This combination of symbol and signal word indicates a video link that is intended to additionally explain a function. In addition, a QR code is displayed to call up the video with a smartphone or tablet.

1.1.1 Always get the latest instructions

If an operating manual or a software manual (hereinafter referred to as 'manual') is changed by the manufacturer, it will be put 'online' immediately. This ensures the compliance of Ecolab Engineering GmbH with the requirements of the 'product monitoring obligation'.

All instructions are provided in  PDF format. To open and view the manuals we recommend the PDF Viewer 'Acrobat' by Adobe (<https://acrobat.adobe.com>).

Through the above measures, Ecolab provides various options for ensuring that you can access the most recent operating instructions at all times

Accessing operating instructions using the website of Ecolab Engineering GmbH

On the manufacturer's website (<https://www.ecolab-engineering.de>), the desired instructions can be searched for and selected under the menu item [Download] / [Operating Instructions.]

Call up instructions with the 'DocuAPP' software for Windows® 10

With the Ecolab 'DocuApp' software for Windows®, all operating instructions, such as catalogues, certificates and CE conformity declarations published by Ecolab Engineering can be downloaded to a Windows® PC (Windows® 10).

To install open the 'Microsoft Store' and enter **DocuAPP** in the search box. The store has the 'DocuApp' software for installation. Follow instructions on your screen for installation.

Accessing operating instructions using a smartphone/tablet

With the Ecolab 'DocuApp' all published operating manuals, catalogues, certificates and CE declarations of conformity from Ecolab Engineering can be accessed with smartphones or tablets (Android & IOS Systems).

The documents shown in the 'DocuApp' are always up-to-date and new versions are displayed immediately. For further information about 'DocuApp', a separate software description ((Art. No. 417102298)) is available.

'Ecolab DocuAPP' instructions downloadable

Download of the software description 'DocuApp' (Article no. 417102298):

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The following describes the installation of 'Ecolab DocuAPP' for 'Android' and 'IOS (Apple)' systems

Installation of 'Ecolab DocuAPP' for Android




On Android phones, the 'Ecolab DocuAPP' can be found in the "Google Play Store".

1. Call up the "Google Play Store" with your Smartphone / Tablet.
2. Enter "Ecolab DocuAPP" in the search box.
3. By entering the search term **Ecolab DocuAPP** together with this symbol, you can find and select the 'Ecolab DocuAPP'.
4. Click on [install].
⇒ 'Ecolab DocuAPP' is being installed.

Via a PC, or browser, the 'Ecolab DocuAPP' can be retrieved by using the following link: <https://play.google.com/store/apps/details?id=ecolab.docuApp>

Installation of 'DocuApp' for IOS (Apple) 

On IOS  smartphones, you can find the '**Ecolab DocuApp**'  in the "APP Store" .

1.  Call up the "APP Store"  with your Smartphone / Tablet.
2.  Go to the search function.
3.  Enter "**Ecolab DocuAPP**" in the search box.
4.  By entering the search term **Ecolab DocuAPP** together with this symbol , you can find and select the '**Ecolab DocuApp**'.
5.  Click on *[install]*.
 ⇒ '**Ecolab DocuApp**'  is being installed.

1.1.2 Symbols, highlights and enumerations

Symbols, safety information

Safety instructions are indicated by symbols in this manual. The safety instructions are preceded by signal words that indicate the degree of hazard.

CAUTION!

This combination of symbol and signal word indicates a potentially dangerous situation that could lead to minor or slight injuries if not avoided.

NOTICE!

This combination of symbol and signal word indicates a potentially dangerous situation that could lead to material damage if not avoided.

ENVIRONMENT!

This combination of symbol and signal word indicates possible dangers to the environment.

The environmental symbol denotes environmental protection measures.

Safety instructions in the operating instructions

Safety instructions can refer to specific, individual operating instructions. Such safety instructions are embedded in the operating instructions, so they do not interrupt the reading flow when executing the action. The signal words described above are used.

Example:

1.  Loosen screw.

2. 

CAUTION!

Risk of trapping on the cover!

Close the cover carefully.

3.  Tighten screw.

Tips and recommendations

This symbol highlights useful tips, recommendations and information for an efficient and trouble-free operation.

Further markings

The following markings are used in this manual to highlight operating instructions, results, collections, references and other elements:

Marking	Explanation
1., 2., 3. ...	Step by step operating instructions
	Results of the operating steps
	References to sections of this manual and related documents
	Collections in no set order
[Button]	Controls (e.g. button, switch), indicators (e.g. signal lights)
'Display'	Screen elements (e.g. buttons, assignment of function keys)

1.1.3 Article numbers / EBS-Article numbers

Both item numbers and EBS numbers could be shown in these operating instructions. EBS numbers are Ecolab-internal item numbers and are used exclusively "internal within the group".

1.1.4 Select copyright

This manual is copyright protected. All rights are reserved by the manufacturer. Making this manual available to third parties, reproduction in any form, even partially, and the exploitation and/or disclosure of the contents without written permission from Ecolab Engineering (hereinafter "the manufacturer") is prohibited except for internal purposes. Any contravention of this will result in claims for damages. The manufacturer reserves the right to assert additional claims.

1.2 Transportation

The measurements of the packaging and the weight of the packaging can be found online Section 12 "Technical data" on page 42.

NOTICE!

Material damage due to improper transportation!

- Transport units can fall or tip over if improperly transported. This can cause a large amount of damage. When unloading the shipping crates, during delivery or even during general shipping, proceed safely and pay attention to the symbols and the information on the packaging.

Danger of putting into operation a piece of transport equipment which has been damaged during transport:

- If damages are discovered during unpacking, do not install or put unit into operation, as otherwise uncontrollable faults can occur.

Transport inspection:

- Examine delivery for completeness and transportation damages and report every damage.
- Damage claims can only be filed within the applicable period for complaints.

In event of transportation damages visible from the outside:

- Do not accept the delivery or accept provisionally.
- Note the extent of damage on transport documents Delivery note of the carrier and initiate a complaint immediately.

1.3 Repairs / returns to Ecolab Engineering GmbH

DANGER!

Conditions for returns

Before returning all parts must be completely free of chemicals inside and outside!

We would like to point out that only clean parts and pumps that have been rinsed with water and are free of dosing agents can be accepted by our customer service in order to eliminate the risk of injury to our personnel from chemical products.

Furthermore, we request that the goods sent in are, if possible due to their size, additionally packed in a suitable bag that prevents residual moisture from escaping into the outer packaging. Enclose a copy of the product data sheet of the dosing chemicals used so that our service staff can prepare for the appropriate use of the PPE.

The return must be requested "online":
<https://www.ecolab-engineering.de/de/kontakt/ruecksendung.html>. Fill in all the details and follow the navigation.

The following documents must be completed:

- Returns form
 - Request the form from Ecolab.
 - Fill out the form correctly and in full.
 - Fill in the declaration of no objection.
 - Send both in advance by fax to: (+49 8662 61-258)
- System components:
 - Free from all impurities (rinsed).
 - Must be dispatched in suitable plastic packaging and in a box in order to avoid any leakage of flushing water.
- Cartons:
 - Addressed to (see):
 - ↳ Chapter 1.8.2 'Technical support contact' on page 10.
 - The word "REPAIR" must be written on a sticker.
 - Include a returns form.

1.4 Packaging

The individual packages are packaged to reflect the expected transport conditions. Only environmentally-friendly materials were used for packaging. The packaging is designed to protect the individual components up to assembly against shipping damage, corrosion and other damage. Do not destroy packaging and only remove it just before assembly.

ENVIRONMENT!


Risk of environmental damage due to incorrect disposal.

Packaging materials are valuable raw materials and can, in many cases, be used again or be usefully processed and recycled.

Incorrect disposal of packaging materials can be a threat to the environment.

- Observe the locally applicable disposal regulations.
- Environmentally-friendly disposal of packaging materials.
- If necessary, hire a specialist to carry out disposal.

Symbols on the packaging

Symbol	Designation	Description
	Top	The sign's arrowheads indicate the top of the package. They must always point upwards, otherwise the contents may be damaged.
	Fragile	Indicates packages with fragile or sensitive content. Handle the package with care, do not drop or knock.
	Keep this product dry	Protect packages from moisture and keep dry.
	Electronic components	Electronic components contained in the package.
	Cold	Protect packages from the cold (frost).
	Stacking	The package may be stacked with other similar packages until the specified maximum number is reached. Pay attention to the exact stacking number.

1.5 Storage

Under certain circumstances, instructions for storage, which go beyond the requirements listed here, can be found on the package. These must be observed accordingly.

Please note the following storage conditions:

- Do not store outdoors.
- Store in a dry and dust-free place.
- Do not expose to aggressive media.
- Protect from sunlight.
- Avoid mechanical vibrations.
- Storage temperature: ± 0 to max. 25 °C.
- Relative humidity: max. 80 %.
- If stored for longer than 3 months, regularly check the general condition of all parts and packaging. If necessary, refresh or renew the preservative.

1.6 Equipment marking – Rating plate

Information on equipment marking or the information on the identification plate can be found in  Section 12 "Technical data" on page 42. Important for all queries is the correct specification of the name and type. This is the only way of ensuring flawless and fast processing.

1.7 Warranty

The manufacturer provides a warranty for operational safety, reliability and performance under the following conditions only:

- Assembly, connection, adjustment, maintenance and repairs must be carried out by qualified and authorised specialists with the aid of the User Manual and all the provided documents.
- The use corresponds to the information in the operating instructions.
- Only original equipment spare parts are used for repairs.

Our products are built, tested and CE certified in accordance with current standards/guidelines. They left the factory in a safe, faultless condition. To reach this condition and to ensure a trouble-free operation, the user must pay attention to all the information / warning information, maintenance regulations, etc. that are contained in all the corresponding Instructions for Use and that may be attached to the product.

The general warranty and service conditions of manufacturer apply.

1.8 Contacts

1.8.1 Manufacturer's contact address

Ecolab Engineering GmbH
Raiffeisenstraße 7
83313 Siegsdorf, Germany

Telephone (+49) 86 62 / 61 0
Fax (+49) 86 62 / 61 166

Email: engineering-mailbox@ecolab.com
<http://www.ecolab-engineering.com>

1.8.2 Technical support contact

ECOLAB Engineering GmbH
Raiffeisenstraße 7
83313 Siegsdorf, Germany
Telephone (+49) 86 62 / 61 234
Fax (+49) 86 62 / 61 166
Email: eursiefb-technicalservice@ecolab.com
<http://www.ecolab-engineering.com>

If you contact technical support, make sure to include the type code in the e-mail. The type code is the only way to identify the metering station and its documentation. You can find the type code on the rating plate.

1.8.3 Returns**Return to:****Ecolab Engineering GmbH
- REPAIR -****Zapfendorfstraße 9
D-83313 Siegsdorf**

Tel.: (+49) 8662 61-0

Fax: (+49) 8662 61-258

Before you send anything back to us, please be sure to read the information under [Chapter 1.3 'Repairs / returns to Ecolab Engineering GmbH'](#) on page 7.

2 Safety

2.1 General safety advice

**WARNING!****Supervision of children handling hand sanitizers**

The licensing for hand sanitizers under German pharmaceutical law is generally granted for grown-ups.

There are no hand sanitizers licensed especially for children.

However, this does not mean that children may not or should not disinfect their hands. Children should however carry out hand disinfecting under supervision. It is important that children do not touch their face with fresh product on their fingers, much less touch their eyes.

When smaller children, it is better if a grown-up disinfects the child's hands by putting the child's hands in between his/her own hands and rubbing on the preparation and waiting for it to absorb (usually 30 seconds).

Alcoholic hand sanitizers contain care substances and are harmless for intact skin. In contrast to soap and water, alcohol changes the natural acidic outside layer of the skin, thereby retaining its skin-protecting properties. The alcohol components of a hand sanitizer containing alcohol evaporate very quickly when applied. Alcohol cannot get into the body via the skin.

2.2 Use for the intended purpose

**NOTICE!**

- The dispensers are suitable for non-contact disinfection of hands, but also for dosing liquid soaps and lotions (for gels please note viscosity).
- For hygienic reasons, the dispenser pumps must be replaced regularly after the (fourth canister change (see ↗ *'Removal of the pump'* on page 32 and ↗ *Chapter 10 'Spare parts and accessories'* on page 38)). The pumps must not be reconditioned, but must be replaced! A spare pump is included in the delivery. The article or EBS number for reordering the pumps can be found in the ↗ *Chapter 10 'Spare parts and accessories'* on page 38.
- The dispensers are designed exclusively for operation with 1.5 V AA Mignon LR6 batteries. Pay attention to the correct polarity (+ / -) of the batteries. Do not use rechargeable batteries.
- Make sure that the pump hose hangs freely in the canister and is not kinked at any point.
- Do not pull on the pump hose, the pump may jam.
- The spray pump is only suitable for thin-bodied, easy-flowing disinfectants without additives (micro/nanoparticles).
- Always pay attention that the device is sufficiently stable.

Used media

- The media that is used should be easily free-flowing at room temperature.
- Check media with higher viscosity (gels, lotions) for suitability.
- Do not use any abrasive or aggressive media.
- Only low-viscosity disinfectants and disinfectant gel (alcoholic or non-alcoholic), liquid soaps and lotions are permitted.
- The system may only be used with products that have been validated by Ecolab.
- The materials/media to be deployed for the intended use of the machine are to be procured and used by the owner/operator.
- The correct handling of these materials/media and the risks involved are the sole responsibility of the Owner.
- Hazard warnings and disposal instructions must be provided by the owner/operator.
- All safety regulations for the handling of chemicals must be maintained and the information contained in the material safety data sheet/product data sheet of the metering medium must be observed.
- Liability is not accepted if invalidated products are used!

Safety data sheets

The safety data sheet is primarily intended for the user so that he or she can take any steps necessary for safeguarding his health and safety at work. Ecolab is well aware of the importance of safety data sheets and the responsibility that they entail. The safety data sheets made available by Ecolab undergo constant control. This is done to ensure that they always contain the latest information.

The best thing to do is to post the safety data sheets right beside the equipment or next to the containers so that the proper countermeasures can be implemented at once in the event of an accident.

Downloading of safety data sheets

Safety data sheets:

The latest respective safety data sheets are available online. To download, go to the link shown below or scan in the illustrated QR code. Once there, you can enter the required product and obtain the corresponding safety data sheet as a download.

<https://safetydata.ecolab.eu/index.php?id=1576&L=1>

Designated sanitizer:

SKINMAN Soft Protect FF	
<p>Anti-viral, especially skin-friendly hand sanitizer for the routine containing Vitamin E, Glycerine & Panthenol</p> <p>Dye and fragrance-free hand sanitizer for all-year use.</p> <p>See also: Appendix A.1 'Ecolab - 116032E-SKINMAN SOFT PROTECT FF' on page 48</p>	

MANODES GP	
<p>Ready-to-use liquid hand sanitizer for the food and beverage industry</p> <p>Wide range of effects, provides optimal protection and long-term safety for several hours, no skin irritation with frequent use, dermatologically tested.</p> <p>See also: ↪ <i>Appendix A.2 'Ecolab - 117783E-MANODES GP' on page 61</i></p>	

Equipment safety

- Protect dispenser against moisture.
- Do not immerse the dispenser in water or expose it to a jet of water (pressure-washer).
- Only fill with the media suitable for the intended use; check potential media for suitability prior to use.
- Do not wash dispenser in the dishwasher.
- Mount drip tray to collect any excess medium so it does not reach the floor.

We accept no liability for damage caused by the following:

- Mechanical damages, especially if the QC seal was damaged
- Damages that can be attributed to improper handling.
- Overvoltage, e.g. Due to using the wrong batteries.
- Using the dispenser for purposes other than the ones approved by us and any damages resulting from that
- Use by persons with limited physical, sensory or mental abilities or a lack of knowledge and experience.

WARNING! Any use that is different to the specified application or any other type of use shall count as incorrect use.

2.2.1 Dangers in handling the device / Important installation instructions

Important warnings for installation Dermados Station and use are printed on the back of the unit, which must be observed to ensure operational safety.

Inflammable substances

Since the disinfectant is alcohol-based, flammability cannot be ruled out.. All measures must therefore be taken to minimise the risk of fire and, if necessary, to fight fires.

Avoid heat exposure, open flame and direct sunlight.

Smoking prohibited

Due to the fire hazard of the sanitizer, there is a general smoking ban exists in the vicinity of the dosing stand in order to prevent ignition.

Explosion protection

To avoid or restrict hazardous explosive atmospheres, it is recommended to provide ventilation permanently or regularly.

Do not use in potentially explosive atmospheres.

Avoid contact with self-flammable substances or strong oxidants.

DANGER!

Risk of slipping

Spilled disinfectant reduces a risk of slipping.

Always wipe the leaked disinfectant immediately and dispose of it properly.

WARNING!

The Dermados Station may tip over during transport or installation!

Pay attention:

- when choosing the place of installation, Dermados Station make sure that it stands on a level surface.
Only on a straight surface is stability guaranteed.
- that it is installed in such a way Dermados Station that it does not narrow or obstruct paths. This prevents tipping over by persons who were "in the way" of the station.
- during transport to the weight and centre of gravity of the Dermados Station. Observe the instructions below and if [Chapter 5.1 'Setup'](#) on page 22 necessary call in a second person for help.

WARNING!

Always wear protective goggles when handling the dosing containers to prevent injury to the eyes through unintentional leakage of disinfectant solution!

WARNING!

In order not to contaminate the interior of the dispenser, we urgently advise you to ensure the necessary sterility and cleanliness in all the steps listed below either by wearing hygienically impeccable gloves (latex) or by disinfecting your hands beforehand.

**CAUTION!**

Make sure that it is placed in a well-ventilated area (10 x / h) Dermados Station so that the inevitable alcoholic vapours can escape.

When using one with canisters, Dermados Station the ventilation slots on the canister cover must never be covered to prevent the risk of fire due to outgassing.

Collect / wipe up possible leaks and spills of the canister with a suitable agent. Do not use automatic systems such as hoovers, suction robots, etc. for this purpose, as the solution containing alcohol can get caught in the appliance and cause it to catch fire.

If a spill is detected in the dosing area or in the storage area, the machine Dermados Station must be shut down until the leak is eliminated or, if necessary, repaired.

**DANGER!**

No mobile phones, heating systems, Dermados Station open lighting systems etc. should be used within 0.5 m of them for explosion protection.

**DANGER!****Risk of burns due to electrostatic charge**

Ethanol burns with a hardly visible flame and can be ignited by static electricity. It is therefore essential to ensure that your hands are dry after using the disinfectant dispenser and before touching another person or object in order to rule out static inflammation.

Ambient conditions

- Only use at relative humidity up to max. 90 % and at room temperature.
- Do not use outdoors, as the permissible temperature ranges (0- 25° C) are easily exceeded or fallen short of.
- No malfunctions due electromagnetic fields (EMC certified).

2.3 Arbitrary conversion and spare parts production**CAUTION!**

Changes or modifications are not permitted without prior, written permission from Ecolab Engineering GmbH and result in the forfeiting of any and all warranty entitlements. Original spare parts and accessories approved by the manufacturer to increase safety. The use of other parts excludes the warranty for the resulting consequences. We would like to point out that the CE conformity expires in case of subsequent conversions

3 Scope of delivery

Presentation	Description	Item no. (EBS no.)
	<p>Dosing stand with dosing dispenser, type Dermados Station (including mounted pump and pre-mounted batteries): For 5 and 20 l canisters (can also be used for 500 ml bottles after conversion)</p> <p>or: For 5 and 20 l canisters (can also be used for 500 ml bottles after conversion), rollable</p> <p>or: For 500-ml bottles</p> <p>or: For 1000-ml bottles</p>	<p>10240138 on request</p> <p>10240173 on request</p> <p>10240166 on request</p> <p>10240167 on request</p>
	Dispenser, single (500 ml)	On request (10012743)
	Dispenser, single (1000 ml)	On request (10012741)
	Key to open dosing dispenser	On request (10015375)
	1 drip tray for the dispenser, 500 ml	On request (10015364)
	1 drip tray for the dispenser, 1000 ml	On request (10015372)
	Spare pump (500 ml, 5 and 20 l) (10 pieces)	On request (10014900)
	or: 1 Replacement pump (1000 ml)	On request (10015196)
	Suction lance for Dermados Station 10240138	On request (on request)
	Stainless steel canister hood	On request (on request)

The appropriate disinfectant containers (bottle or canister) are supplied with the respective dosing stands so that the units are "ready to use" (see Chapter 10.1 'Dosing container' on page 39 and Appendix A 'Safety data sheets' on page 48).

4 Functional description

This is a Dermados Station stand-alone device with touch-free operation for safe application of hand disinfectants.

The dispensers are also suitable for dosing liquid soaps and lotions. The integrated sensor detects an approaching hand and triggers an automatic dispensing process.

It Dermados Station is used in entrance areas of office buildings, hotels, hospitals, food and beverage manufacturers, cosmetics, pharmaceutical and other industrial companies as well as schools, kindergartens and supermarkets.


It is a battery-powered system that can be set up autonomously.

It Dermados Station is available in different versions, so that it can be operated with a 5 or 20 l canister or with 500 and 1000 ml bottles, depending on the version (see also [Chapter 4.1 'Variants' on page 19](#)).

The 5 and 20 l versions are particularly suitable for highly frequented disinfection processes, as the large volume of the disinfection container means that canisters do not need to be replaced as often. With the 20 l version, over 6000 dosages are possible without canister or battery change.

As a further version of the 5 or 20 l version, a "rollable" version is optionally available. This version facilitates installation and is particularly advantageous if the station is to be relocated.

General set-up

- 1 Housing corpus**
 - a Lock
 - b ON / OFF switch (+ special function)
 - c Adjustment button for dispensing quantity (+ special function)
 - d Dispenser pump, writable
 - e Battery compartment
 - f Sensor mounting plate
 - g LED
 - h Nozzle
 - i Dosing hose to nozzle (short)
 - j Dosing hose to canister (long)
 - k Cylinder holder
- 2 Housing cover**
 - l Latch
- 3 Drip tray**
- 4 Floor-mounted console**
- 5 Dosing agent box for 5 and 20 l canisters**

4.1 Variants

Dosing stand with dosing dispenser (500 ml, 5 litres and 20 l)

Presentation	Description	Item no. (EBS no.)
 <p>500 ml 5,000 ml / 20,000 ml</p> <p>500 ml / 5,000 ml / 20,000 ml rollbar / scrollable / enroulable</p>	<p>Complete dosing stand with dosing dispenser Including mounted pump and pre-mounted batteries. Also available as rollable version. Suitable disinfectant containers: 500 ml bottle, 5 and 20 l canister</p>	<p>10240138 on request</p> <p>1024173 on request</p>

The use of a 500 ml bottle requires a change-over including pump change (🔗 'Conversion from 5 or 20 l canister variant to 500 ml bottle variant' on page 41).

Dosing stand with dosing dispenser (500 ml)

Presentation	Description	Item no. (EBS no.)
	<p>Complete dosing stand with dosing dispenser Including mounted pump and pre-mounted batteries. Suitable disinfectant containers: 500 ml Flasche</p>	<p>10240166 on request</p>

Dosing stand with dosing dispenser (1000 ml)

Presentation	Description	Item no. (EBS no.)
	<p>Complete dosing stand with dosing dispenser Including mounted pump and pre-mounted batteries. Also available as rollable version. Suitable disinfectant containers: 1000 ml cylinder</p>	<p>10240167 on request</p> <p>10240140 on request</p>

Dispenser (500 ml)

Presentation	Description	Item no. (EBS no.)
	<p>Dispenser (500 ml) incl. mounted pump and pre-mounted batteries.</p> <p>Suitable disinfectant containers: <i>500 ml cylinder</i></p>	<p>10240151 on request</p>

Dispenser (1000 ml)

Presentation	Description	Item no. (EBS no.)
	<p>Dispenser (1000 ml) incl. mounted pump and pre-mounted batteries.</p> <p>Suitable disinfectant containers: <i>1000 ml cylinder</i></p>	<p>10240140 on request</p>

5 Installation and assembly

- Personnel: ■ Trained person
 ■ Service personnel
- Protective equipment: ■ Protective eyewear
 ■ Chemical-resistant protective gloves

WARNING!

In order not to contaminate the interior of the dispenser, we urgently advise you to ensure the necessary sterility and cleanliness in all the steps listed below either by wearing hygienically impeccable gloves (latex) or by disinfecting your hands beforehand.

DANGER!

Risk of burns due to electrostatic charge

Ethanol burns with a hardly visible flame and can be ignited by static electricity. It is therefore essential to ensure that your hands are dry after using the disinfectant dispenser and before touching another person or object in order to rule out static inflammation.

DANGER!

No mobile phones, heating systems, Dermados Station open lighting systems etc. should be used within 0.5 m of them for explosion protection.

CAUTION!

Make sure that it is placed in a well-ventilated area (10 x / h) Dermados Station so that the inevitable alcoholic vapours can escape.

When using one with canisters, Dermados Station the ventilation slots on the canister cover must never be covered to prevent the risk of fire due to outgassing.

Collect / wipe up possible leaks and spills of the canister with a suitable agent. Do not use automatic systems such as hoovers, suction robots, etc. for this purpose, as the solution containing alcohol can get caught in the appliance and cause it to catch fire.

If a spill is detected in the dosing area or in the storage area, the machine Dermados Station must be shut down until the leak is eliminated or, if necessary, repaired.

5.1 Setup

WARNING!

The Dermados Station may tip over during transport or installation!

Pay attention:

- when choosing the place of installation, Dermados Station make sure that it stands on a level surface.
Only on a straight surface is stability guaranteed.
- that it is installed in such a way Dermados Station that it does not narrow or obstruct paths. This prevents tipping over by persons who were "in the way" of the station.
- during transport to the weight and centre of gravity of the Dermados Station. Observe the instructions below and if Chapter 5.1 'Setup' on page 22 necessary call in a second person for help.

- Set up the ward in a well-ventilated area.
- Place it in a clearly visible position so that it can be recognised and used.
- Make sure that no paths are blocked to prevent accidents.
- Always keep them Dermados Station ready for operation.
It is the operator's responsibility to ensure proper functioning.
- We strongly recommend to store both brand name batteries (4 x 1,5 V AA Mignon LR6), as well as dosing agents.

Roller version

In order to simplify installation, there is a "roller version" of the Dermados Station which makes transport very easy, especially when used with 5 and 20 l canisters (see Table on page 17).

5.2 Assembly

Variant: 5 and 20 l canisters


Fig. 1: Dermados Station Canister assembly / change


- | | |
|---------------------|---------------------------------|
| 1 Long suction hose | 6 Stainless steel canister hood |
| 2 Suction pipe | 7 canister |
| 3 Foot valve | 8 Suction lance canister lid |
| 4 Base panel | 9 Safety instructions |
| 5 Lock device | |

Mount the dosing container (canister) and connect it to the suction lance

1.  **CAUTION!**
Mount the dosing container (canister) and Dermados Station connect it to the suction lance (pos. 9)

Remove the pre-assembled stainless steel canister bonnet (pos. 6) upwards.

- 2.** Place the canister in the centre of the base plate (pos. 4)
- 3.** Remove the lid of the canister.
The suction lance (pos. 2) is fitted with a matching canister lid (pos. 8).
- 4.** Insert the suction lance into the canister (pos. 7) as far as possible.
- 5.** Put on the suction lance canister lid (pos. 8) and tighten it hand-tight.

 **DANGER!**
All safety regulations for handling chemicals must always be observed (see also  Appendix A 'Safety data sheets' on page 48)!

- 6.** Slide the stainless steel canister cover onto the brackets from above.

7.  *We recommend to lock the canister cover with a suitable lock (Pos. 5) from unauthorised access. The rollable version is secured with the enclosed lateral Allen screws.*

Dermados Station: 500 and 1000 ml cylinders

The variants for using 500 and 1000 ml bottles differ in the size of the dispenser. The 500 ml bottle can therefore not be used in the dispenser of the 1000 ml version (article no. 10240167), nor can the 1000 ml in the 500 ml version (article no. 10240166).

Insert disinfectant bottle

1. ▶ Open dosing dispenser with key (🔑 Chapter 3 'Scope of delivery' on page 17).
2. ▶ Remove the cap from the 500 or 1000 ml bottle (pos. 4).
3. ▶ Insert filled bottle. To do this, guide the bottle over the suction hose from below and ensure that the hose hangs centrally in the bottle and is not kinked.
4. ▶ Place the bottle neck over the guide spigot (on the underside of the sensor carrier plate).
5. ▶ Place the bottle on the bottle stop (pos. 5).
6. ▶ Remove the key.
7. ▶ Close dosing dispenser.

6 Commissioning

- Personnel:
- Trained person
 - Service personnel
- Protective equipment:
- Protective eyewear
 - Chemical-resistant protective gloves

WARNING!

In order not to contaminate the interior of the dispenser, we urgently advise you to ensure the necessary sterility and cleanliness in all the steps listed below either by wearing hygienically impeccable gloves (latex) or by disinfecting your hands beforehand.

DANGER!

No mobile phones, heating systems, Dermados Station open lighting systems etc. should be used within 0.5 m of them for explosion protection.

DANGER!

Risk of burns due to electrostatic charge

Ethanol burns with a hardly visible flame and can be ignited by static electricity. It is therefore essential to ensure that your hands are dry after using the disinfectant dispenser and before touching another person or object in order to rule out static inflammation.

Preparations

1. ▶ Open dispenser with key (see [Chapter 3 'Scope of delivery'](#) on page 17).

Fig. 2: Peristaltic pump

2. ▶ Label the pump with the installation date (Pos. 1).

Fig. 3: Insert batteries (4 x AA Mignon LR6)

3. ▶ The batteries are already inserted by the factory.
If this is not the case, proceed as follows:
 - Open battery compartment (Pos. 2).
 - Insert four new 1.5 V AA Mignon LR6 batteries (item 3), making sure the polarity is correct (see imprint in the battery compartment lid).
 - Close battery compartment (Pos. 2).
4. ▶ Remove the key.
5. ▶ Close dispenser.

Preparing for dosing

1. ▶ Hold one hand (or both hands) in the sensor area.
The sensitive range is factory set at the centre of the dispenser.

For setting see: [Chapter 3 'Set dosing level \(dosing quantity\):'](#) on page 27

2. ▶ Remove hand and hold in sensor area once again.
 - ⇒ The dispenser doses again.
 - ⇒ After two or three doses, the pump is filled and the dispenser is ready for use.

Set dosing level (dosing quantity):

Set dosing level (dosing quantity):

https://www.youtube.com/watch?v=oP_XXo85Zfg

The discharge quantity per dosing process can be adjusted in four steps with the housing cover open by pressing the round setting knob, pos. 1 (approx. One second).

Factory setting: Step 2 = the dispenser doses about 3.2 ml of liquid.

The entire preset dispensing quantity is dosed without interruptions. From Level 2 onwards, we recommend using both hands to catch the dosing amount.

The dosing amount can be less if the media has high viscosity.


Fig. 4: Setting the dosing level (dosing amount)

1 Adjusting knob 'Dosing quantity'

2 LED

1. ➤ The unit is preset to Level 2 by the factory. The dispenser doses approx. 3.2 ml liquid. Press the setting knob (position 1) to activate step 3.
 - ⇒ The red LED (pos. 2) flashes three times = the dispenser doses approx. 4.8 ml of liquid.
2. ➤ Press the setting knob again to activate level 4.
 - ⇒ The red LED flashes four times = the dispenser doses approx. 6.4 ml liquid.

From here the settings level starts again from the beginning: From level 1 to level 4.

3. ➤ Press the button again to return to dosing level 1
 - ⇒ The red LED flashes once = the dispenser doses approx. 1.6 ml of liquid.

7 Operation - Hand sanitizing

WARNING!

Supervision of children handling hand sanitizers

The licensing for hand sanitizers under German pharmaceutical law is generally granted for grown-ups.

There are no hand sanitizers licensed especially for children.

However, this does not mean that children may not or should not disinfect their hands. Children should however carry out hand disinfecting under supervision. It is important that children do not touch their face with fresh product on their fingers, much less touch their eyes.

When smaller children, it is better if a grown-up disinfects the child's hands by putting the child's hands in between his/her own hands and rubbing on the preparation and waiting for it to absorb (usually 30 seconds).

Alcoholic hand sanitizers contain care substances and are harmless for intact skin. In contrast to soap and water, alcohol changes the natural acidic outside layer of the skin, thereby retaining its skin-protecting properties. The alcohol components of a hand sanitizer containing alcohol evaporate very quickly when applied. Alcohol cannot get into the body via the skin.

DANGER!

No mobile phones, heating systems, Dermados Station open lighting systems etc. should be used within 0.5 m of them for explosion protection.

DANGER!

Risk of burns due to electrostatic charge

Ethanol burns with a hardly visible flame and can be ignited by static electricity. It is therefore essential to ensure that your hands are dry after using the disinfectant dispenser and before touching another person or object in order to rule out static inflammation.

As long as the hand is held in the sensor area, the system doses once. If the dosage is to be triggered several times, the hand must be removed and then brought back into the sensor area. If higher dosing quantities are required, it is easier and more sensible to increase the dosing levels. To change the dosing level, see ☞ 'Set dosing level (dosing quantity):' on page 27.

Make sure that the dispenser is ready for operation.

An LED is mounted on the top of dispenser, easily visible:

- **LED lights up blue** = Dispenser is ready for use.
- **Red LED is visible** = malfunction, e.g. Replace batteries
- **LED lights up blue-red** (mixed colour) = electronic malfunction or incorrect operation. Perform a "reset" as a first measure (see ☞ Chapter 8 'Description and rectification of malfunctions' on page 30). If this does not remedy the problem, please ask the manufacturer (☞ Chapter 1.8 'Contacts' on page 10) whether you can remedy the fault yourself.

Hand sanitizing

Fig. 5: Hand sanitizing

1. ▶ Hold one or both hands in the sensor area.
 ⇨ The preset dosing amount is dispensed.

2. ▶

NOTICE!

Rub the disinfectant for at least 30 seconds and be sure to wet all areas of the hands.

DANGER!

Risk of burns due to electrostatic charge

Ethanol burns with a hardly visible flame and can be ignited by static electricity. It is therefore essential to ensure that your hands are dry after using the disinfectant dispenser and before touching another person or object in order to rule out static inflammation.

8 Description and rectification of malfunctions

The dosing dispenser is subject to regulatory requirements (applicable to area of EU: Maschinen- Richtlinie, ElektroG, BatterieG).

Therefore you should not attempt doing any repairs yourself, but contact your dealer. You may possibly be able to rectify some breakdowns.

Fault description	Cause	Remedy
LED alternately blinks red and blue.	Square and round button was pressed accidentally (e.g. while inserting the batteries).	Perform a reset. To do this, remove a battery, then press down on the square (ON OFF) button and hold it down until the removed battery is reinserted.
	Defective electronics.	Please contact vendor.
Red LED is blinking (after inserting new batteries).	Batteries are not correctly inserted / twisted.	Check position of batteries, carefully check polarity.
Dispenser is not dosing.	Canister is empty.	Insert new canister (see ↗ <i>'Exchange canister of dosing agent'</i> on page 34).
	Pump hose is kinked.	Check hose and rectify kink, if possible.
	Pump hose is sticking to the canister wall (when bottle is almost empty).	Examine and position hose in centre of canister, possible replace canister (see ↗ <i>'Exchange canister of dosing agent'</i> on page 34).
	Wrong medium (high viscosity).	Examine medium for suitability.
Dispenser dosing quantity is too much.	Use more fluid medium.	Set to highest dosing level (see ↗ <i>'Set dosing level (dosing quantity):'</i> on page 27).
	Set to dosing level 1 = approx. 1.6 ml.	Set to dosing level 2 = approx. 3.2 ml (see ↗ <i>'Set dosing level (dosing quantity):'</i> on page 27).
Pump runs sluggishly.	Pump in use for more than 6 months.	Replace pump (possibly every 3-4 months) (see ↗ <i>'Removal of the pump'</i> on page 32).
	Use more fluid medium.	Examine medium for suitability.
Sensor is not working, although hand is being held under the sensor.	Hand held down too low.	Bring hand closer to sensor.
	Triggering distance needs to be adjusted.	Recalibrate sensor (see ↗ <i>'Sensor calibration'</i> on page 37).
	Defective sensor/ electronics.	Please contact vendor.
Difficulty completely closing dispenser.	Key is still in lock.	Remove key.
	Hose is kinked.	Install hose accurately.

9 Service

- Personnel: ■ Trained person
 ■ Service personnel
- Protective equipment: ■ Protective eyewear
 ■ Chemical-resistant protective gloves

WARNING!

In order not to contaminate the interior of the dispenser, we urgently advise you to ensure the necessary sterility and cleanliness in all the steps listed below either by wearing hygienically impeccable gloves (latex) or by disinfecting your hands beforehand.

DANGER!

Risk of burns due to electrostatic charge

Ethanol burns with a hardly visible flame and can be ignited by static electricity. It is therefore essential to ensure that your hands are dry after using the disinfectant dispenser and before touching another person or object in order to rule out static inflammation.

DANGER!

No mobile phones, heating systems, Dermados Station open lighting systems etc. should be used within 0.5 m of them for explosion protection.

- *Every fourth canister change (but at least every **three** months) should be cleaned with a damp cloth dipped in warm soapy water or the dispenser housing disinfected:*
 - *Do not immerse the appliance in water or wash it with a high-pressure cleaner, for example.*
 - *Dispenser cannot be accidentally triggered if the front of the housing is opened prior to cleaning.*
 - *Do not use any caustic (chlorinated) or abrasive cleaners.*
- *If the pump hose is blocked after standing around for a longer period of time, then using a new pump is the easier and cheaper solution when compared with cleaning.*
- *The device's guarantee lapses if the QC seal is damaged!*

9.1 Service

Pump exchange:

<https://youtu.be/IKHx5F-hnwY>

CAUTION!

Dispenser care

Inadequate dispenser hygiene can result in microbial colonisation of the dispenser system. The most frequent germs are aerobic sporulating organisms, *Corynebacteria*, various cocci and *pseudomonas*. Hand sanitizer and soap residues in the hose or on the nozzle can also cause problems.

NOTICE!

For hygienic reasons, the dispenser pumps are designed as disposable pumps which must be replaced at regular intervals every 3 - 4 months, or after 4 canister changes at the latest.

The pumps must not be reconditioned, but must be replaced!

Label pump with date prior to using it for first time!

When used for the first time or after a long period of use, the tube may be stuck in the pump head at the point where the rollers pinch it off. This can be remedied by moving the visible part of the gear on the back of the pump head back and forth several times.

Removal of the pump


Fig. 6: Remove the pump

1. ▶ Open dispenser with key (see [Chapter 3 'Scope of delivery'](#) on page 17).
2. ▶ Take the nozzle out of the bracket (Pos. 1).
3. ▶ Turn the pump (pos. 3) clockwise (approx. 20°) until the lock is released.
4. ▶ Remove the pump.

Installation of the pump

i A suitable pump is always included in the delivery (see ↪ Chapter 10 'Spare parts and accessories' on page 38). If this has already been used, we strongly recommend that you keep at least **one pump in stock so that you always have a spare pump available** for upcoming maintenance work. You can order suitable pumps as spare parts (↪ Chapter 10 'Spare parts and accessories' on page 38)


Fig. 7: Re-install the pump

1. ➤ Label the new pump (Pos. 1) with the current date, insert it and then turn it counter-clockwise (Pos. 2), until it engages.
2. ➤ **500 or 1000 ml variant:**
Insert the longer end of the silicone hose (pos. 3), which comes out of the pump (suction side), into the bottle, latch the shorter end (pressure side) with the nozzle attached to the tip of the sensor carrier plate (pos. 4).
5 or 20 l variant:
Put the longer end of the silicone hose (pos. 3) coming from the pump (suction side) over the hose connector and connect it to the hose leading into the canister. Latch the shorter end (print side) with the nozzle attached to the tip of the sensor carrier plate (pos. 4).

NOTICE!

Note on the hose connector

Please make sure that the hose ends are pulled at least 12 mm over the hose connector to avoid loosening of the hose connection (e.g. due to transport vibrations or pump movements)!

Hose length of the spare pump for 5 and 20 l version

Please note that when ordering replacement pumps, the length of the suction-side pump hose may have to be adjusted before it is pushed over the hose connector.

3. ➤ If properly positioned, the nozzle can be reinserted without any effort required.
4. ➤ Remove the key and close the dispenser.
5. ➤ Operate the dispenser until the used medium emerges from the nozzle.

If this does not work, please use the instructions at ↪ Chapter 8 'Description and rectification of malfunctions' on page 30.

⇒ Dispenser is ready for use!

Exchange canister of dosing agent

Replace dosing agent canister:

<https://youtu.be/p7COXWkyHag>

NOTICE!

To ensure an uninterrupted supply to the system, you should replace the canister in good time and, if necessary, keep others in reserve.

CAUTION!

Whenever canisters are changed, the person responsible must ensure that the cap is correctly screwed on and that any leaked liquid is immediately removed professionally.

1. ▶ Open and remove the security lock, if present, on the canister cover.
2. ▶ Lift canister lid off and outwards out of the console.
3. ▶ Pull the suction lance out of the canister.
Make sure that no residual dosing agent drips onto the floor or wipe it up again immediately with a suitable cloth.
4. ▶ The dosing container (canister) is supplied with a pre-drilled lid suitable for the suction lance. This lid is to be replaced with the original sealing lid of the dosing container.
5. ▶ Place the new canister in the middle of the base plate.
6. ▶ Insert the suction lance through the hole in the canister lid into the dosing container until the foot valve touches the container base.
7. ▶ Slide the canister cover onto the brackets from above.
8. ▶ We recommend protecting the canister cover from unauthorised access with a suitable lock.
9. ▶ Perform two or three pump operations to refill the hose system.
When the pump dispenser is open, no unintended triggering of a pump operation is possible (the ON/OFF switch is deactivated).

Replace disinfectant bottle


Replace the disinfectant bottle:
In Work process

NOTICE!

To ensure uninterrupted supply to the system, you should replace the 500 or 1000 ml disinfectant bottles in good time and keep them in reserve if necessary.

Bottle change

1. ➤ Open dosing dispenser with key (↪ *Chapter 3 'Scope of delivery' on page 17*).
2. ➤ Pull the dosing hose (pos. 1), coming from the pump, out of the bottle.
3. ➤ Remove the cap of the new bottle (pos. 4).
4. ➤ Guide the new bottle over the suction hose from below and ensure that the hose hangs centrally in the bottle and is not kinked.
5. ➤ Place the bottle neck over the guide spigot (on the underside of the sensor carrier plate).
6. ➤ Lift the bottle stopper (pos. 5) so that the bottle is clamped.
7. ➤ Remove the key.
8. ➤ Close dosing dispenser.

Battery replacement

Battery replacement

<https://youtu.be/YPf51vVM5us>

The operation of the dispensers requires very little power. This is a set of **4 x 1.5 V AA Mignon LR6** brand-name batteries suffices for about two years (approx. 30,000 dosing cycles). As a rule, the batteries last even longer.

NOTICE!

- Never mix new and used batteries.
- Always use batteries of the same type.
- In the event the dispenser has not been used for a longer time, the batteries should be removed to prevent potential loss of charge or leakage.
- Do not use rechargeable batteries, they do not reach the required voltage even when fully charged.

Inserting batteries

1. ▶ Open dosing dispenser with key (☞ *Chapter 3 'Scope of delivery' on page 17*).
2. ▶ Depending on the stand variant, either remove the bottle or detach the suction-side pump hose from the hose connector. Pass the pump hose next to each pump housing to access the battery compartment.

Fig. 8: Insert batteries (4 x AA Mignon LR6)

3. ▶ Open battery compartment (Pos. 2).
4. ▶ Insert 4 new 1.5 V AA Mignon LR6 batteries (Pos. 3), pay attention to the proper polarity during the process (see labelling on the battery compartment lid).
5. ▶ Close the battery compartment cover (pos. 2).
6. ▶ Either replace the bottle (500 or 1000 ml) or connect the pump hose to the canister hose (5 or 20 l variant) via the hose connector.
7. ▶ Close dispenser housing.

Sensor calibration

Sensor calibration
In Work process

1. ➤ Open (🔑 *Chapter 3 'Scope of delivery' on page 17*) the dispenser housing with the key and open the lid.
2. ➤ Depending on the stand variant, either remove the bottle or detach the suction-side pump hose from the hose connector. Pass the pump hose next to each pump housing to access the battery compartment.
3. ➤ Open the battery compartment and remove one of the four batteries or, for initial installation, insert three batteries.
4. ➤ Using two fingers of one hand (e.g. the index finger and middle finger), press the ON/OFF switch (square; pos. 1) and the delivery rate setting knob (round; pos. 2) simultaneously.
5. ➤ Use your other hand to install the fourth battery, keeping both buttons pressed down.
6. ➤ The LED (pos. 3) lights up 1x blue and then 3 x blue again in rapid succession.
7. ➤ Now release the two buttons. Close the dispenser housing quickly.
8. ➤ Place your hand or an object (e.g. white cardboard) in the sensor beam (item 4) at the position to which the trigger distance is to be reset. As soon as the sensor has registered the hand or object, the LED flashes in two colours, slows down and finally flashes once in blue.
⇒ The setting procedure is completed!
9. ➤ Remove the hand/object from the sensor beam.
10. ➤ Open the dispenser housing again and close the battery compartment.
11. ➤ Either replace the bottle (500 or 1000 ml bottle variant) or connect the pump hose to the canister hose (5 or 20 l variant) via the hose connector.
12. ➤ Check that the desired release distance is set.

10 Spare parts and accessories

Figure	Name	Article No.	EBS. No.
	Peristaltic squeeze pump Dermados TF 500	10240152	On request
	Peristaltic squeeze pump Dermados TF 1000	10240145	On request
	Dermados TouchLess Dispenser 500 ml	10240151	On request
	Dermados TouchLess Dispenser 1000 ml	10240140	On request
	Dermados Sticker Skinman Soft	10240154	On request
	P-Dermados Sticker Manodes GP	10240135	On request
	Straight hose connector, PVDF, 4 mm	415100681	10032889
	Suction lance set Consisting of: <ul style="list-style-type: none"> ■ Foot valve (non-return valve), for a longer service life of the peristaltic pump ■ Screw connection on canister lid ■ complete hose (from pump to foot valve) 	10240200	On request
	Foot filter PP	10240144	On request
	Hose, \varnothing 3 x 1.5 (3/6) EVA natural	417400911	On request

10.1 Dosing container

Figure	Name	Article No.	EBS. No.
Skinman Soft Protect FF: Virucidal liquid hand disinfectant with vitamin E, glycerine and panthenol			
	500 ml	On request	On request
	1000 ml	On request	On request
	5 l	On request	On request
	20 l	On request	On request
Manodes GP: Liquid hand disinfectant for all areas except patient treatment and care in healthcare.			
	500 ml	On request	On request
	1000 ml	On request	On request
	5 l	On request	On request
	20 l	On request	On request

Order of hand disinfectants via:
Ordermanagement@ecolab.com

11 Conversion - canister to bottle or reactivation of canister version

- Personnel:
- Trained person
 - Service personnel
- Protective equipment:
- Protective eyewear
 - Chemical-resistant protective gloves

WARNING!

In order not to contaminate the interior of the dispenser, we urgently advise you to ensure the necessary sterility and cleanliness in all the steps listed below either by wearing hygienically impeccable gloves (latex) or by disinfecting your hands beforehand.

DANGER!

Risk of burns due to electrostatic charge

Ethanol burns with a hardly visible flame and can be ignited by static electricity. It is therefore essential to ensure that your hands are dry after using the disinfectant dispenser and before touching another person or object in order to rule out static inflammation.

DANGER!

No mobile phones, heating systems, Dermados Station open lighting systems etc. should be used within 0.5 m of them for explosion protection.

A conversion is only possible if the existing disinfection station is equipped with a 500 ml dispenser. It is therefore important to know which type of dispenser is installed on your console so that the desired container fits.

Conversion from 5 or 20 l canister variant to 500 ml bottle variant

When converting the 5 and 20 l canister variant to a 500 ml variant, minor adjustments must be made.

Modify hose routing

1. ➤ Open dosing dispenser with key (↪ *Chapter 3 'Scope of delivery' on page 17*).
⇒ The hose guide in the 20 l sanitizer tank becomes visible (see figure on left).
2. ➤ Remove the short dosing hose (pos. 1), coming from the pump, (Pos. 2) from the hose connector.
3. ➤ Remove the long suction line including hose connector (pos. 3) from the dispenser housing.
4. ➤ Remove the cap from the 500 ml bottle (pos. 4).
5. ➤ Guide the new bottle over the suction hose from below and ensure that the hose hangs centrally in the bottle and is not kinked.
6. ➤ Place the bottle neck over the guide spigot (on the underside of the sensor carrier plate).
7. ➤ Lift the bottle stop (pos. 5) so that the bottle is firmly clamped.
8. ➤ Remove the key.
9. ➤ Close dosing dispenser.

Conversion from 500 ml bottle variant to 5 or 20 l canister variant

*To reactivate the dispenser with the 5 or 20 l disinfection container, reconnect the short dosing hose (pos. 1) to the hose connector (pos. 2) and make the connection ↪ *Chapter 5.2 'Assembly' on page 23* as described below.*

12 Technical specifications

In this chapter you will find the technical data of the Dermados Station.

Specification	Value
Housing material incl. mounting plate	ABS-plastic, impact-resistant
Pump material	ABS plastic
Hose material	Silicon
Colour: Housing and mounting plate	white/grey
Batteries	4 x 1,5 V AA Mignon LR6

Data	Value	Unit
Voltage	6	V
Protection rating	IP 22	
Weight of console (single)	approx. 19	kg
Weight of dosing dispenser (single and excl. batteries)	approx. 0.85	kg
Total weight	approx. 20	kg
Deformation temperature	> 85	°C

12.1 Type plate

The rating plate is located at the top on the reverse side of the Dermados station.


Fig. 9: Type plate

- | | |
|------------------------------------|--------------------------|
| 1 Type designation | 4 Production code |
| 2 Serial number | 5 CE conformity |
| 3 Specification of the noise level | 6 Manufacturer's address |

12.2 Dimensions

Data	Value	Unit
Dimensions including packaging	1350 x 385 x 68	mm
Measurements excluding packaging	1312 x 295 x 372	mm

Dispenser (single)


Fig. 10: Dimensions dispenser, single

Dermados Station 500 ml and 1000 ml


Fig. 11: Dimensions Dermados Station 500 ml and 1000 ml

Dermados Station 5 and 20 I


Fig. 12: Dimensions Dermados Station 5 and 20 I

13 Decommissioning, disassembly, environmental protection

Withdrawal from service

1. ➤ Open dosing dispenser with key (☞ *Chapter 3 'Scope of delivery' on page 17*).
2. ➤ Open the battery compartment.
3. ➤ Remove batteries.

Disassembly

DANGER!

- Keep the workplace tidy and clean. Loosely stacked components and tools (or leaving components and tools lying around) can cause accidents.
- Disassemble the components professionally.

1. ➤ Remove dosing agent canister or bottle and close with original cap.
2. ➤ Open dosing dispenser with key (☞ *Chapter 3 'Scope of delivery' on page 17*).
3. ➤ Pull the dispenser upwards out of the holder.
4. ➤ Unscrew the holding plate of the dosing dispenser from the stand.
5. ➤ Dermados Station into its individual parts.

Disposal and environmental protection

ENVIRONMENT!

Risk of environmental damage due to incorrect disposal!

Do not dispose of any components in the domestic waste. Take them instead to the appropriate collection points for recycling.

All components are to be disposed of in accordance with prevailing local environmental regulations. Dispose of them accordingly, depending on the nature of the product, existing regulations and with due regard for current provisions and criteria. We ask you to observe Directive 2012/19/EU on waste electrical and electronic devices, that has the objective and purpose of reducing and/or avoiding the waste from recyclable raw materials. This directive requires member states of the EU to increase the collection rate of electronic waste so that it can be recycled.

If in doubt, contact your local authority, or an approved waste disposal service provider, for information on correct disposal.

Recycle the dismantled components:

- Scrap all metals.
- Electrical waste and electronic components must be recycled.
- Recycle all plastic parts.
- Dispose of all other components in line with their material characteristics.
- Hand in batteries at communal collection points or dispose of them through a specialist.

14 EC Declaration/Declaration of Conformity

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Annex 1 to WI-EU-RDE-602 Rev. 0 / 2019-06-13

Fig. 13: Conformity according to EC Machinery Directive 2006/42/EC | EMC Directive 2004/108/EC.

i Due to technical changes, it is possible that the 'Declaration of Conformity / CE declaration' may change. The most recent 'Declaration of conformity / CE declaration' will therefore be published on the Internet.

Appendix

A Safety data sheets**A.1 Ecolab - 116032E-SKINMAN SOFT PROTECT FF**

Name	Data
Designation	Safety Data Sheet 116032E-SKINMAN SOFT PROTECT FF
Type	SKINMAN SOFT PROTECT FF
Number	
Type of manual	Safety data sheet
Manufacturer	Ecolab Deutschland GmbH Ecolab-Allee 1 40789 Monheim am Rhein

Skinman Soft Protect FF**Section: 1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING****1.1 Product identifier**

Product name : Skinman Soft Protect FF

Product code : 116032E

Use of the
Substance/Mixture : Hand Sanitizer

Substance type: : Mixture

For professional users only.

Product dilution information : No dilution information provided.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Skin disinfectant

Recommended restrictions : Reserved for industrial and professional use.
on use**1.3 Details of the supplier of the safety data sheet**Company : Ecolab Deutschland GmbH
Ecolab-Allee 1
40789 Monheim am Rhein, Germany +49 (0)2173 599 0
OfficeService.DEDUS@ecolab.com**1.4 Emergency telephone number**Emergency telephone : +4932221096286
number : +32-(0)3-575-5555 Trans-EuropeanPoison Information Centre : +49 (0)551 38318854
telephone number

Date of Compilation/Revision : 06.04.2020

Version : 3.2

Section: 2. HAZARDS IDENTIFICATION**2.1 Classification of the substance or mixture****Classification (REGULATION (EC) No 1272/2008)**Flammable liquids, Category 2 H225
Chronic aquatic toxicity, Category 3 H412**2.2 Label elements****Labelling (REGULATION (EC) No 1272/2008)**

Skinman Soft Protect FF

Hazard pictograms	:		
Signal Word	:	Danger	
Hazard Statements	:	H225 H412	Highly flammable liquid and vapour. Harmful to aquatic life with long lasting effects.
Precautionary Statements	:	Prevention: P210 P273	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid release to the environment.

2.3 Other hazards

None known.

Section: 3. COMPOSITION/INFORMATION ON INGREDIENTS

3.2 Mixtures

Hazardous components

Chemical Name	CAS-No. EC-No. REACH No.	Classification REGULATION (EC) No 1272/2008	Concentration: [%]
ethanol	64-17-5 200-578-6 01-2119457610-43	Flammable liquids Category 2; H225	>= 50 - <= 100
Myristyl Alcohol	112-72-1 204-000-3 01-2119485910-33	Eye irritation Category 2; H319 Chronic aquatic toxicity Category 1; H410	>= 1 - < 2.5
Substances with a workplace exposure limit :			
butanone	78-93-3 201-159-0 01-2119457290-43	Flammable liquids Category 2; H225 Eye irritation Category 2; H319 Specific target organ toxicity - single exposure Category 3; H336	>= 0.5 - < 1
glycerin	56-81-5 200-289-5 01-2119471987-18	Not Classified;	>= 0.25 - < 0.5

For the full text of the H-Statements mentioned in this Section, see Section 16.

Section: 4. FIRST AID MEASURES

4.1 Description of first aid measures

In case of eye contact	:	Rinse with plenty of water.
In case of skin contact	:	Rinse with plenty of water.
If swallowed	:	Rinse mouth. Get medical attention if symptoms occur.
If inhaled	:	Get medical attention if symptoms occur.

Skinman Soft Protect FF

4.2 Most important symptoms and effects, both acute and delayed

See Section 11 for more detailed information on health effects and symptoms.

4.3 Indication of immediate medical attention and special treatment needed

Treatment : Treat symptomatically.

Section: 5. FIREFIGHTING MEASURES

5.1 Extinguishing media

Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media : High volume water jet

5.2 Special hazards arising from the substance or mixture

Specific hazards during firefighting : Fire Hazard
Keep away from heat and sources of ignition.
Flash back possible over considerable distance.
Beware of vapours accumulating to form explosive concentrations.
Vapours can accumulate in low areas.

Hazardous combustion products : Depending on combustion properties, decomposition products may include following materials:
Carbon oxides

5.3 Advice for firefighters

Special protective equipment for firefighters : Use personal protective equipment.

Further information : Use water spray to cool unopened containers. Fire residues and contaminated fire extinguishing water must be disposed of in accordance with local regulations. In the event of fire and/or explosion do not breathe fumes.

Section: 6. ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel : Remove all sources of ignition. Ensure clean-up is conducted by trained personnel only. Refer to protective measures listed in sections 7 and 8.

Advice for emergency responders : If specialised clothing is required to deal with the spillage, take note of any information in Section 8 on suitable and unsuitable materials.

6.2 Environmental precautions

Environmental precautions : Do not allow contact with soil, surface or ground water.

6.3 Methods and materials for containment and cleaning up

Skinman Soft Protect FF

Methods for cleaning up : Eliminate all ignition sources if safe to do so. Stop leak if safe to do so. Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13). Flush away traces with water. For large spills, dike spilled material or otherwise contain material to ensure runoff does not reach a waterway.

6.4 Reference to other sections

See Section 1 for emergency contact information.
 For personal protection see section 8.
 See Section 13 for additional waste treatment information.

Section: 7. HANDLING AND STORAGE

7.1 Precautions for safe handling

Advice on safe handling : Use only with adequate ventilation. Keep away from fire, sparks and heated surfaces. Take necessary action to avoid static electricity discharge (which might cause ignition of organic vapours). Wash hands thoroughly after handling. Open drum carefully as content may be under pressure. In case of mechanical malfunction, or if in contact with unknown dilution of product, wear full Personal Protective Equipment (PPE).

Hygiene measures : Handle in accordance with good industrial hygiene and safety practice. Remove and wash contaminated clothing before re-use. Wash face, hands and any exposed skin thoroughly after handling.

7.2 Conditions for safe storage, including any incompatibilities

Requirements for storage areas and containers : Keep away from heat and sources of ignition. Keep in a cool, well-ventilated place. Keep away from oxidizing agents. Keep out of reach of children. Keep container tightly closed. Store in suitable labeled containers.

Storage temperature : 0 °C to 25 °C

7.3 Specific end uses

Specific use(s) : Skin disinfectant

Section: 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

Occupational Exposure Limits

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
ethanol	64-17-5	AGW	200 ppm 380 mg/m3	TRGS 900
Further information	DFG	Senate commission for the review of compounds at the work place dangerous for the health (MAK-commission).		
	Y	When there is compliance with the OEL and biological tolerance values, there is no risk of harming the unborn child		
Myristyl Alcohol	112-72-1	AGW (Vapour and	20 ppm	TRGS 900

Skinman Soft Protect FF

		aerosols)	178 mg/m3	
Further information	AGS	Commission for dangerous substances		
	11	Sum of vapor and aerosols.		
butanone	78-93-3	AGW	200 ppm 600 mg/m3	TRGS 900
Further information	DFG	Senate commission for the review of compounds at the work place dangerous for the health (MAK-commission).		
	EU	European Union (The EU has established a limit value: deviations in value and peak limit are possible)		
	H	Skin absorption		
	Y	When there is compliance with the OEL and biological tolerance values, there is no risk of harming the unborn child		
glycerin	56-81-5	AGW (Inhalable fraction)	200 mg/m3	TRGS 900
Further information	DFG	Senate commission for the review of compounds at the work place dangerous for the health (MAK-commission).		
	Y	When there is compliance with the OEL and biological tolerance values, there is no risk of harming the unborn child		

Biological occupational exposure limits

Substance name	CAS-No.	Control parameters	Sampling time	Basis
Alkyl ketone	Proprietary Ingredient	2-butanone: 2 mg/l (Urine)	Immediately after exposition or after working hours	TRGS 903

8.2 Exposure controls

Appropriate engineering controls

Engineering measures : Effective exhaust ventilation system. Maintain air concentrations below occupational exposure standards.

Individual protection measures

Hygiene measures : Handle in accordance with good industrial hygiene and safety practice. Remove and wash contaminated clothing before re-use. Wash face, hands and any exposed skin thoroughly after handling.

Eye/face protection (EN 166) : No special protective equipment required.

Hand protection (EN 374) : No special protective equipment required.

Skin and body protection (EN 14605) : No special protective equipment required.

Respiratory protection (EN 143, 14387) : None required if airborne concentrations are maintained below the exposure limit listed in Exposure Limit Information. Use certified respiratory protection equipment meeting EU requirements(89/656/EEC, (EU) 2016/425), or equivalent, when respiratory risks cannot be avoided or sufficiently limited by technical means of collective protection or by measures, methods or procedures of work organization.

Environmental exposure controls

General advice : Consider the provision of containment around storage vessels.

Section: 9. PHYSICAL AND CHEMICAL PROPERTIES

Skinman Soft Protect FF

9.1 Information on basic physical and chemical properties

Appearance	: liquid
Colour	: clear, colourless
Odour	: alcohol-like
pH	: 6.0 - 7.5, 100 %
Flash point	: 16 °C closed cup
Odour Threshold	: Not applicable and/or not determined for the mixture
Melting point/freezing point	: Not applicable and/or not determined for the mixture
Initial boiling point and boiling range	: Not applicable and/or not determined for the mixture
Evaporation rate	: Not applicable and/or not determined for the mixture
Flammability (solid, gas)	: Not applicable and/or not determined for the mixture
Upper explosion limit	: Not applicable and/or not determined for the mixture
Lower explosion limit	: Not applicable and/or not determined for the mixture
Vapour pressure	: Not applicable and/or not determined for the mixture
Relative vapour density	: Not applicable and/or not determined for the mixture
Relative density	: 0.81 - 0.83
Water solubility	: soluble
Solubility in other solvents	: Not applicable and/or not determined for the mixture
Partition coefficient: n-octanol/water	: Not applicable and/or not determined for the mixture
Auto-ignition temperature	: Not applicable and/or not determined for the mixture
Thermal decomposition	: Not applicable and/or not determined for the mixture
Viscosity, kinematic	: Not applicable and/or not determined for the mixture
Explosive properties	: Not applicable and/or not determined for the mixture
Oxidizing properties	: The substance or mixture is not classified as oxidizing.

9.2 Other information

Not applicable and/or not determined for the mixture

Section: 10. STABILITY AND REACTIVITY

10.1 Reactivity

No dangerous reaction known under conditions of normal use.

10.2 Chemical stability

Stable under normal conditions.

10.3 Possibility of hazardous reactions

No dangerous reaction known under conditions of normal use.

10.4 Conditions to avoid

Skinman Soft Protect FF

Heat, flames and sparks.

10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products

Depending on combustion properties, decomposition products may include following materials:
Carbon oxides

Section: 11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects

Information on likely routes of exposure : Inhalation, Eye contact, Skin contact

Product

Acute oral toxicity : There is no data available for this product.

Acute inhalation toxicity : There is no data available for this product.

Acute dermal toxicity : There is no data available for this product.

Skin corrosion/irritation : There is no data available for this product.

Serious eye damage/eye irritation : There is no data available for this product.

Respiratory or skin sensitization : There is no data available for this product.

Carcinogenicity : There is no data available for this product.

Reproductive effects : There is no data available for this product.

Germ cell mutagenicity : There is no data available for this product.

Teratogenicity : There is no data available for this product.

STOT - single exposure : There is no data available for this product.

STOT - repeated exposure : There is no data available for this product.

Aspiration toxicity : There is no data available for this product.

Components

Acute oral toxicity : ethanol
LD50 rat: 10,470 mg/kg

butanone
LC50 rat: 2,193 mg/kg
Test substance: Information given is based on data obtained from similar substances.

Skinman Soft Protect FF

glycerin
LD50 rat: 18,300 mg/kg

Components

Acute inhalation toxicity : ethanol
4 h LC50 rat: 117 mg/l
Test atmosphere: vapour

butanone
4 h LC50 rat: 34.4 mg/l
Test atmosphere: vapour

Components

Acute dermal toxicity : ethanol
LD50 rabbit: > 15,800 mg/kg

Myristyl Alcohol
LD50 rabbit: 8,000 mg/kg

butanone
LD50 rat: > 8,050 mg/kg

glycerin
LD50 rabbit: 23,000 mg/kg

Potential Health Effects

Eyes : Health injuries are not known or expected under normal use.

Skin : Health injuries are not known or expected under normal use.

Ingestion : Health injuries are not known or expected under normal use.

Inhalation : Health injuries are not known or expected under normal use.

Chronic Exposure : Health injuries are not known or expected under normal use.

Health injuries are not known or expected under normal use.

Experience with human exposure

Eye contact : No symptoms known or expected.

Skin contact : No symptoms known or expected.

Ingestion : No symptoms known or expected.

Inhalation : No symptoms known or expected.

Section: 12. ECOLOGICAL INFORMATION

12.1 Ecotoxicity

Environmental Effects : Harmful to aquatic life with long lasting effects.

Skinman Soft Protect FF

Product

Toxicity to fish : no data available

Toxicity to daphnia and other aquatic invertebrates : no data available

Toxicity to algae : no data available

Components

Toxicity to fish : ethanol
96 h LC50 Pimephales promelas (fathead minnow): > 100 mg/l

Myristyl Alcohol
96 h LC50 Oncorhynchus mykiss (rainbow trout): > 1 mg/l

butanone
96 h LC50 Pimephales promelas (fathead minnow): 2,993 mg/l

glycerin
96 h LC50 Fish: 855 mg/l

Components

Toxicity to daphnia and other aquatic invertebrates : Myristyl Alcohol
48 h EC50 Daphnia magna (Water flea): 3.2 mg/l

butanone
48 h EC50 Daphnia magna (Water flea): 308 mg/l

Components

Toxicity to algae : Myristyl Alcohol
96 h EC50 Desmodesmus subspicatus (green algae): > 10 mg/l

butanone
96 h EC50 Pseudokirchneriella subcapitata (algae): 2,029 mg/l

12.2 Persistence and degradability

Product

no data available

Components

Biodegradability : ethanol
Result: Readily biodegradable.

Myristyl Alcohol
Result: Readily biodegradable.

butanone
Result: Readily biodegradable.

glycerin
Result: Readily biodegradable.

12.3 Bioaccumulative potential

no data available

Skinman Soft Protect FF

12.4 Mobility in soil

no data available

12.5 Results of PBT and vPvB assessment

Product

Assessment : This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

12.6 Other adverse effects

no data available

Section: 13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with the European Directives on waste and hazardous waste. Waste codes should be assigned by the user, preferably in discussion with the waste disposal authorities.

13.1 Waste treatment methods

- Product : The product should not be allowed to enter drains, water courses or the soil. Where possible recycling is preferred to disposal or incineration. If recycling is not practicable, dispose of in compliance with local regulations. Dispose of wastes in an approved waste disposal facility.
- Contaminated packaging : Dispose of as unused product. Empty containers should be taken to an approved waste handling site for recycling or disposal. Do not re-use empty containers. Dispose of in accordance with local, state, and federal regulations.
- Guidance for Waste Code selection : Organic wastes containing dangerous substances. If this product is used in any further processes, the final user must redefine and assign the most appropriate European Waste Catalogue Code. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable European (EU Directive 2008/98/EC) and local regulations.

Section: 14. TRANSPORT INFORMATION

The shipper/consignor/sender is responsible to ensure that the packaging, labeling, and markings are in compliance with the selected mode of transport.

Land transport (ADR/ADN/RID)

- 14.1 UN number : 1170
- 14.2 UN proper shipping name : ETHANOL SOLUTION
- 14.3 Transport hazard class(es) : 3
- 14.4 Packing group : II
- 14.5 Environmental hazards : No

Skinman Soft Protect FF

14.6 Special precautions for user : None

Air transport (IATA)

14.1 UN number : 1170
14.2 UN proper shipping name : Ethanol solution
14.3 Transport hazard class(es) : 3
14.4 Packing group : II
14.5 Environmental hazards : No
14.6 Special precautions for user : None

Sea transport (IMDG/IMO)

14.1 UN number : 1170
14.2 UN proper shipping name : ETHANOL SOLUTION
14.3 Transport hazard class(es) : 3
14.4 Packing group : II
14.5 Environmental hazards : No
14.6 Special precautions for user : None
14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code : Not applicable.

Section: 15. REGULATORY INFORMATION

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

National Regulations

Take note of Dir 94/33/EC on the protection of young people at work.

Hazard class for water : WGK 1
Classification according to AwSV, Annex 1

German storage class : 3

15.2 Chemical Safety Assessment

No Chemical Safety Assessment has been carried out on the product.

Section: 16. OTHER INFORMATION

Procedure used to derive the classification according to REGULATION (EC) No 1272/2008

Classification	Justification
Flammable liquids 2, H225	Based on product data or assessment
Chronic aquatic toxicity 3, H412	Calculation method

Full text of H-Statements

H225 Highly flammable liquid and vapour.
H319 Causes serious eye irritation.

Skinman Soft Protect FF

H336 May cause drowsiness or dizziness.
H410 Very toxic to aquatic life with long lasting effects.

Full text of other abbreviations

ADN – European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR – European Agreement concerning the International Carriage of Dangerous Goods by Road; AICS – Australian Inventory of Chemical Substances; ASTM – American Society for the Testing of Materials; bw – Body weight; CLP – Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008; CMR – Carcinogen, Mutagen or Reproductive Toxicant; DIN – Standard of the German Institute for Standardisation; DSL – Domestic Substances List (Canada); ECHA – European Chemicals Agency; EC-Number – European Community number; ECx – Concentration associated with x% response; ELx – Loading rate associated with x% response; EmS – Emergency Schedule; ENCS – Existing and New Chemical Substances (Japan); ErCx – Concentration associated with x% growth rate response; GHS – Globally Harmonized System; GLP – Good Laboratory Practice; IARC – International Agency for Research on Cancer; IATA – International Air Transport Association; IBC – International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 – Half maximal inhibitory concentration; ICAO – International Civil Aviation Organization; IECSC – Inventory of Existing Chemical Substances in China; IMDG – International Maritime Dangerous Goods; IMO – International Maritime Organization; ISHL – Industrial Safety and Health Law (Japan); ISO – International Organisation for Standardization; KECI – Korea Existing Chemicals Inventory; LC50 – Lethal Concentration to 50 % of a test population; LD50 – Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL – International Convention for the Prevention of Pollution from Ships; n.o.s. – Not Otherwise Specified; NO(A)EC – No Observed (Adverse) Effect Concentration; NO(A)EL – No Observed (Adverse) Effect Level; NOELR – No Observable Effect Loading Rate; NZIoC – New Zealand Inventory of Chemicals; OECD – Organization for Economic Co-operation and Development; OPPTS – Office of Chemical Safety and Pollution Prevention; PBT – Persistent, Bioaccumulative and Toxic substance; PICCS – Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR – (Quantitative) Structure Activity Relationship; REACH – Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID – Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT – Self-Accelerating Decomposition Temperature; SDS – Safety Data Sheet; TCSI – Taiwan Chemical Substance Inventory; TRGS – Technical Rule for Hazardous Substances; TSCA – Toxic Substances Control Act (United States); UN – United Nations; vPvB – Very Persistent and Very Bioaccumulative

Prepared by : Regulatory Affairs

Numbers quoted in the MSDS are given in the format: 1,000,000 = 1 million and 1,000 = 1 thousand. 0.1 = 1 tenth and 0.001 = 1 thousandth

REVISED INFORMATION: Significant changes to regulatory or health information for this revision is indicated by a bar in the left-hand margin of the SDS.

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

A.2 Ecolab - 117783E-MANODES GP

Name	Data
Designation	Safety data sheet 117783E-MANODES GP
Type	MANODES GP
Number	
Type of manual	Safety data sheet
Manufacturer	Ecolab Deutschland GmbH Ecolab-Allee 1 40789 Monheim am Rhein

MANODES GP**Section: 1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING****1.1 Product identifier**

Product name : MANODES GP
Product code : 117783E
Use of the Substance/Mixture : Hand Sanitizer
Substance type: : Mixture

For professional users only.

Product dilution information : No dilution information provided.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Professional Hand Cleaner /skin disinfectants
Recommended restrictions on use : Reserved for industrial and professional use.

1.3 Details of the supplier of the safety data sheet

Company : Ecolab Deutschland GmbH
Ecolab-Allee 1
40789 Monheim am Rhein, Germany +49 (0)2173 599 0
OfficeService.DEDUS@ecolab.com

1.4 Emergency telephone number

Emergency telephone number : +4932221096286
+32-(0)3-575-5555 Trans-European
Poison Information Centre telephone number : +49 (0)551 38318854

Date of Compilation/Revision : 14.04.2020
Version : 1.1

Section: 2. HAZARDS IDENTIFICATION**2.1 Classification of the substance or mixture****Classification (REGULATION (EC) No 1272/2008)**

Flammable liquids, Category 2 H225

2.2 Label elements

Labelling (REGULATION (EC) No 1272/2008)

MANODES GP

Hazard pictograms	:		
Signal Word	:	Danger	
Hazard Statements	:	H225	Highly flammable liquid and vapour.
Precautionary Statements	:	Prevention: P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

2.3 Other hazards

None known.

Section: 3. COMPOSITION/INFORMATION ON INGREDIENTS

3.2 Mixtures

Hazardous components

Chemical Name	CAS-No. EC-No. REACH No.	Classification REGULATION (EC) No 1272/2008	Concentration: [%]
ethanol	64-17-5 200-578-6 01-2119457610-43	Flammable liquids Category 2; H225	>= 50 - <= 100
Substances with a workplace exposure limit :			
glycerin	56-81-5 200-289-5 01-2119471987-18	Not Classified;	>= 1 - < 2.5
butanone	78-93-3 201-159-0 01-2119457290-43	Flammable liquids Category 2; H225 Eye irritation Category 2; H319 Specific target organ toxicity - single exposure Category 3; H336	>= 0.5 - < 1

For the full text of the H-Statements mentioned in this Section, see Section 16.

Section: 4. FIRST AID MEASURES

4.1 Description of first aid measures

In case of eye contact	:	Rinse with plenty of water.
In case of skin contact	:	Rinse with plenty of water.
If swallowed	:	Rinse mouth. Get medical attention if symptoms occur.
If inhaled	:	Get medical attention if symptoms occur.

4.2 Most important symptoms and effects, both acute and delayed

See Section 11 for more detailed information on health effects and symptoms.

4.3 Indication of immediate medical attention and special treatment needed

MANODES GP

Treatment : Treat symptomatically.

Section: 5. FIREFIGHTING MEASURES

5.1 Extinguishing media

Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media : High volume water jet

5.2 Special hazards arising from the substance or mixture

Specific hazards during firefighting : Fire Hazard
Keep away from heat and sources of ignition.
Flash back possible over considerable distance.
Beware of vapours accumulating to form explosive concentrations.
Vapours can accumulate in low areas.

Hazardous combustion products : Depending on combustion properties, decomposition products may include following materials:
Carbon oxides

5.3 Advice for firefighters

Special protective equipment for firefighters : Use personal protective equipment.

Further information : Use water spray to cool unopened containers. Fire residues and contaminated fire extinguishing water must be disposed of in accordance with local regulations. In the event of fire and/or explosion do not breathe fumes.

Section: 6. ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel : Remove all sources of ignition. Refer to protective measures listed in sections 7 and 8.

Advice for emergency responders : If specialised clothing is required to deal with the spillage, take note of any information in Section 8 on suitable and unsuitable materials.

6.2 Environmental precautions

Environmental precautions : Do not allow contact with soil, surface or ground water.

6.3 Methods and materials for containment and cleaning up

Methods for cleaning up : Eliminate all ignition sources if safe to do so. Stop leak if safe to do so. Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13). Flush away traces with water. For large spills, dike spilled material or otherwise contain material to ensure runoff does not reach a waterway.

MANODES GP

6.4 Reference to other sections

See Section 1 for emergency contact information.
 For personal protection see section 8.
 See Section 13 for additional waste treatment information.

Section: 7. HANDLING AND STORAGE

7.1 Precautions for safe handling

Advice on safe handling : Keep away from fire, sparks and heated surfaces. Take necessary action to avoid static electricity discharge (which might cause ignition of organic vapours). Open drum carefully as content may be under pressure.

7.2 Conditions for safe storage, including any incompatibilities

Requirements for storage areas and containers : Keep away from heat and sources of ignition. Keep in a cool, well-ventilated place. Keep away from oxidizing agents. Keep out of reach of children. Keep container tightly closed. Store in suitable labeled containers.

Storage temperature : 0 °C to 25 °C

7.3 Specific end uses

Specific use(s) : Professional Hand Cleaner /skin desinfectants

Section: 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

Occupational Exposure Limits

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
ethanol	64-17-5	AGW	200 ppm 380 mg/m3	TRGS 900
Further information	DFG	Senate commission for the review of compounds at the work place dangerous for the health (MAK-commission).		
	Y	When there is compliance with the OEL and biological tolerance values, there is no risk of harming the unborn child		
glycerin	56-81-5	AGW (Inhalable fraction)	200 mg/m3	TRGS 900
Further information	DFG	Senate commission for the review of compounds at the work place dangerous for the health (MAK-commission).		
	Y	When there is compliance with the OEL and biological tolerance values, there is no risk of harming the unborn child		
butanone	78-93-3	AGW	200 ppm 600 mg/m3	TRGS 900
Further information	DFG	Senate commission for the review of compounds at the work place dangerous for the health (MAK-commission).		
	EU	European Union (The EU has established a limit value: deviations in value and peak limit are possible)		
	H	Skin absorption		
	Y	When there is compliance with the OEL and biological tolerance values, there is no risk of harming the unborn child		

Biological occupational exposure limits

MANODES GP

Substance name	CAS-No.	Control parameters	Sampling time	Basis
Alkyl ketone	Proprietary Ingredient	2-butanone: 2 mg/l (Urine)	Immediately after exposition or after working hours	TRGS 903

DNEL

Hydrogen peroxide	:	End Use: Workers Exposure routes: Inhalation Potential health effects: Short-term - local Value: 3 mg/m ³
	:	End Use: Workers Exposure routes: Inhalation Potential health effects: Long-term local effects Value: 1.4 mg/m ³

8.2 Exposure controls

Appropriate engineering controls

Engineering measures : Good general ventilation should be sufficient to control worker exposure to airborne contaminants.

Individual protection measures

Eye/face protection (EN 166) : No special protective equipment required.

Hand protection (EN 374) : No special protective equipment required.

Skin and body protection (EN 14605) : No special protective equipment required.

Respiratory protection (EN 143, 14387) : None required if airborne concentrations are maintained below the exposure limit listed in Exposure Limit Information. Use certified respiratory protection equipment meeting EU requirements(89/656/EEC, (EU) 2016/425), or equivalent, when respiratory risks cannot be avoided or sufficiently limited by technical means of collective protection or by measures, methods or procedures of work organization.

Environmental exposure controls

General advice : Consider the provision of containment around storage vessels.

Section: 9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

Appearance : liquid
 Colour : colourless
 Odour : alcohol-like
 pH : 5
 Flash point : 17 °C
 Odour Threshold : Not applicable and/or not determined for the mixture

MANODES GP

Melting point/freezing point	: Not applicable and/or not determined for the mixture
Initial boiling point and boiling range	: > 35 °C
Evaporation rate	: Not applicable and/or not determined for the mixture
Flammability (solid, gas)	: Not applicable and/or not determined for the mixture
Upper explosion limit	: Not applicable and/or not determined for the mixture
Lower explosion limit	: Not applicable and/or not determined for the mixture
Vapour pressure	: Not applicable and/or not determined for the mixture
Relative vapour density	: Not applicable and/or not determined for the mixture
Relative density	: 0.84 - 0.87
Water solubility	: soluble
Solubility in other solvents	: Not applicable and/or not determined for the mixture
Partition coefficient: n-octanol/water	: Not applicable and/or not determined for the mixture
Auto-ignition temperature	: Not applicable and/or not determined for the mixture
Thermal decomposition	: Not applicable and/or not determined for the mixture
Viscosity, kinematic	: Not applicable and/or not determined for the mixture
Explosive properties	: Not applicable and/or not determined for the mixture
Oxidizing properties	: Not applicable and/or not determined for the mixture

9.2 Other information

Not applicable and/or not determined for the mixture

Section: 10. STABILITY AND REACTIVITY

10.1 Reactivity

No dangerous reaction known under conditions of normal use.

10.2 Chemical stability

Stable under normal conditions.

10.3 Possibility of hazardous reactions

No dangerous reaction known under conditions of normal use.

10.4 Conditions to avoid

Heat, flames and sparks.

10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products

Depending on combustion properties, decomposition products may include following materials:
Carbon oxides

MANODES GP

Section: 11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects

Information on likely routes of exposure : Inhalation, Eye contact, Skin contact

Product

Acute oral toxicity : There is no data available for this product.

Acute inhalation toxicity : There is no data available for this product.

Acute dermal toxicity : There is no data available for this product.

Skin corrosion/irritation : There is no data available for this product.

Serious eye damage/eye irritation : There is no data available for this product.

Respiratory or skin sensitization : There is no data available for this product.

Carcinogenicity : There is no data available for this product.

Reproductive effects : There is no data available for this product.

Germ cell mutagenicity : There is no data available for this product.

Teratogenicity : There is no data available for this product.

STOT - single exposure : There is no data available for this product.

STOT - repeated exposure : There is no data available for this product.

Aspiration toxicity : There is no data available for this product.

Components

Acute oral toxicity : ethanol
LD50 rat: 10,470 mg/kg

glycerin
LD50 rat: 18,300 mg/kg

butanone
LC50 rat: 2,193 mg/kg

Test substance: Information given is based on data obtained from similar substances.

Components

Acute inhalation toxicity : ethanol
4 h LC50 rat: 117 mg/l
Test atmosphere: vapour

MANODES GP

butanone
4 h LC50 rat: 34.4 mg/l
Test atmosphere: vapour

Components

Acute dermal toxicity : ethanol
LD50 rabbit: > 15,800 mg/kg

glycerin
LD50 rabbit: 23,000 mg/kg

butanone
LD50 rat: > 8,050 mg/kg

Potential Health Effects

Eyes : Health injuries are not known or expected under normal use.

Skin : Health injuries are not known or expected under normal use.

Ingestion : Health injuries are not known or expected under normal use.

Inhalation : Health injuries are not known or expected under normal use.

Chronic Exposure : Health injuries are not known or expected under normal use.

Experience with human exposure

Eye contact : No symptoms known or expected.

Skin contact : No symptoms known or expected.

Ingestion : No symptoms known or expected.

Inhalation : No symptoms known or expected.

Section: 12. ECOLOGICAL INFORMATION

12.1 Ecotoxicity

Environmental Effects : This product has no known ecotoxicological effects.

Product

Toxicity to fish : no data available

Toxicity to daphnia and other aquatic invertebrates : no data available

Toxicity to algae : no data available

Components

Toxicity to fish : ethanol
96 h LC50 Pimephales promelas (fathead minnow): > 100 mg/l

glycerin
96 h LC50 Fish: 855 mg/l

butanone

MANODES GP

96 h LC50 Pimephales promelas (fathead minnow): 2,993 mg/l

Components

Toxicity to daphnia and other aquatic invertebrates : butanone
48 h EC50 Daphnia magna (Water flea): 308 mg/l

Components

Toxicity to algae : butanone
96 h EC50 Pseudokirchneriella subcapitata (algae): 2,029 mg/l

12.2 Persistence and degradability

Product

no data available

Components

Biodegradability : ethanol
Result: Readily biodegradable.

glycerin
Result: Readily biodegradable.

butanone
Result: Readily biodegradable.

12.3 Bioaccumulative potential

no data available

12.4 Mobility in soil

no data available

12.5 Results of PBT and vPvB assessment

Product

Assessment : This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

12.6 Other adverse effects

no data available

Section: 13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with the European Directives on waste and hazardous waste. Waste codes should be assigned by the user, preferably in discussion with the waste disposal authorities.

13.1 Waste treatment methods

Product : Where possible recycling is preferred to disposal or incineration. If recycling is not practicable, dispose of in compliance with local regulations. Dispose of wastes in an approved waste disposal

MANODES GP

facility.

Contaminated packaging : Dispose of as unused product. Empty containers should be taken to an approved waste handling site for recycling or disposal. Do not re-use empty containers. Dispose of in accordance with local, state, and federal regulations.

Guidance for Waste Code selection : Organic wastes containing dangerous substances. If this product is used in any further processes, the final user must redefine and assign the most appropriate European Waste Catalogue Code. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable European (EU Directive 2008/98/EC) and local regulations.

Section: 14. TRANSPORT INFORMATION

The shipper/consignor/sender is responsible to ensure that the packaging, labeling, and markings are in compliance with the selected mode of transport.

Land transport (ADR/ADN/RID)

14.1 UN number : 1170
14.2 UN proper shipping name : ETHANOL SOLUTION
14.3 Transport hazard class(es) : 3
14.4 Packing group : II
14.5 Environmental hazards : No
14.6 Special precautions for user : None

Air transport (IATA)

14.1 UN number : 1170
14.2 UN proper shipping name : Ethanol solution
14.3 Transport hazard class(es) : 3
14.4 Packing group : II
14.5 Environmental hazards : No
14.6 Special precautions for user : None

Sea transport (IMDG/IMO)

14.1 UN number : 1170
14.2 UN proper shipping name : ETHANOL SOLUTION
14.3 Transport hazard class(es) : 3
14.4 Packing group : II
14.5 Environmental hazards : No
14.6 Special precautions for user : None
14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code : Not applicable.

MANODES GP**Section: 15. REGULATORY INFORMATION****15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture****National Regulations**

Take note of Dir 94/33/EC on the protection of young people at work.

Hazard class for water : WGK 1
Classification according to AwSV, Annex 1

German storage class : 3

15.2 Chemical Safety Assessment

No Chemical Safety Assessment has been carried out on the product.

Section: 16. OTHER INFORMATION**Procedure used to derive the classification according to REGULATION (EC) No 1272/2008**

Classification	Justification
Flammable liquids 2, H225	Based on product data or assessment

Full text of H-Statements

H225 Highly flammable liquid and vapour.
H319 Causes serious eye irritation.
H336 May cause drowsiness or dizziness.

Full text of other abbreviations

ADN – European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR – European Agreement concerning the International Carriage of Dangerous Goods by Road; AICS – Australian Inventory of Chemical Substances; ASTM – American Society for the Testing of Materials; bw – Body weight; CLP – Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008; CMR – Carcinogen, Mutagen or Reproductive Toxicant; DIN – Standard of the German Institute for Standardisation; DSL – Domestic Substances List (Canada); ECHA – European Chemicals Agency; EC-Number – European Community number; ECx – Concentration associated with x% response; ELx – Loading rate associated with x% response; EmS – Emergency Schedule; ENCS – Existing and New Chemical Substances (Japan); ErCx – Concentration associated with x% growth rate response; GHS – Globally Harmonized System; GLP – Good Laboratory Practice; IARC – International Agency for Research on Cancer; IATA – International Air Transport Association; IBC – International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 – Half maximal inhibitory concentration; ICAO – International Civil Aviation Organization; IECSC – Inventory of Existing Chemical Substances in China; IMDG – International Maritime Dangerous Goods; IMO – International Maritime Organization; ISHL – Industrial Safety and Health Law (Japan); ISO – International Organisation for Standardization; KECI – Korea Existing Chemicals Inventory; LC50 – Lethal Concentration to 50 % of a test population; LD50 – Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL – International Convention for the Prevention of Pollution from Ships; n.o.s. – Not Otherwise Specified; NO(A)EC – No Observed (Adverse) Effect Concentration; NO(A)EL – No Observed (Adverse) Effect Level; NOELR – No Observable Effect Loading Rate; NZIoC – New Zealand Inventory of Chemicals; OECD – Organization for Economic Co-operation and Development; OPPTS – Office of Chemical Safety and Pollution Prevention; PBT – Persistent, Bioaccumulative and Toxic substance; PICCS – Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR – (Quantitative) Structure Activity Relationship; REACH – Regulation (EC) No 1907/2006 of the European

MANODES GP

Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID – Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT – Self-Accelerating Decomposition Temperature; SDS – Safety Data Sheet; TCSI – Taiwan Chemical Substance Inventory; TRGS – Technical Rule for Hazardous Substances; TSCA – Toxic Substances Control Act (United States); UN – United Nations; vPvB – Very Persistent and Very Bioaccumulative

Prepared by : Regulatory Affairs

Numbers quoted in the MSDS are given in the format: 1,000,000 = 1 million and 1,000 = 1 thousand. 0.1 = 1 tenth and 0.001 = 1 thousandth

REVISED INFORMATION: Significant changes to regulatory or health information for this revision is indicated by a bar in the left-hand margin of the SDS.

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

Dokumenten-Nr.:	Dermados Station
document no.:	
Erstelldatum:	23.09.2020
date of issue:	
Version / Revision:	MAN048961, Rev. 2-09.2020
version / revision:	
Letze Änderung:	16.09.2020
last changing:	

Copyright [Ecolab Engineering GmbH](#), 2020

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Reproduction, also in part, only with permission of
[Ecolab Engineering GmbH](#)

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Alle Rechte vorbehalten *All rights reserved*
der Firma [Ecolab Engineering GmbH](#)