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## Editorial

Dear readers,

summer is slowly coming to an end and the time of cozy autumn evenings to relax and lose yourself in a book begins. Why not add aseptica to your holiday reading list?

In this issue, you can read all about the latest developments in the worlds of hygiene and reprocessing and find out everything that you, as a professional, need to know about these and other topics.

In our Latest News section, Dr Kaufmann discusses the potential for widespread testing to become a key element of our fight against the pandemic, including a closer look at the pros and cons of the various test types and their limitations.

In his article on “Requirements for the validation of cleaning and disinfection processes – DIN 58341”, Mr Kruse covers the new standard and the opportunities and benefits it offers for validators and operators.

In the Technology and Hygiene section, you can learn more about the interaction between reprocessing products and reprocessing equipment, with a focus on lumened dental instruments and the associated adapters for washer-disinfectors.

I hope you enjoy reading this issue of aseptica.

Stay healthy!

Stella Nehr-Werner

## Contents

### Latest News

Coronavirus SARS-CoV-2: Testing as a key element of the fight against the pandemic 27

New Instrument Reprocessing Working Group brochure in the series: “Instrument Reprocessing: Reprocessing of Instruments to Retain Value” 31

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“Miele MOVE” portal monitors laundry machines 39

## Report

### Fitness trackers detect long COVID

The fitness trackers that many people now wear day and night can be used to work out the average duration of diseases. In a US study published in JAMA Network Open (2021; <https://jamanetwork.com/journals/jamanetworkopen/fullarticle/2781687>; 10.1001/jamanetworkopen.2021.15959), researchers found that heart rates remained elevated in patients diagnosed with COVID-19 for a period of two to three months, while sleep and physical activity returned to normal levels faster. The study also found new evidence of relative bradycardia – a slight drop in heart rate – in the acute phase of the disease.

The 37,146 US participants in the study submitted their data between 25 March 2020 and 24 January 2021 to the “<https://www.scripps.edu/science-and-medicine/translational-institute/>” Research Translational Institute, where the data was evaluated. 827 of the participants told the app that they had been tested for SARS-CoV-2 due to acute respiratory symptoms. The test was positive for 234 participants.

Source: aertzeblatt.de

[www.aseptica.com](http://www.aseptica.com)  
Download a digital copy of the latest edition now and browse through the extensive archive.

### Technology & Hygiene

Requirements for the validation of cleaning and disinfection processes – DIN 58341 40

The interaction between reprocessing products (specifically lumened dental instruments) and reprocessing equipment (adapters in washer-disinfectors) 43

### Miscellaneous & Legal Notice

“Three questions for ...” Dr Hubert Holz 46

New advisory body member: Ines Konschake 46

Event announcements 47  
25. German Society for Sterile Supplies & 21st WFHSS Congress

# Coronavirus SARS-CoV-2: Testing as a key element of the fight against the pandemic

Sabine Kaufmann

In her best-selling book “Pale Rider: The Spanish Flu of 1918 and How It Changed the World” (published by Jonathan Cape), Laura Spinney paints a multifaceted picture of how this pandemic changed society. Over 100 years ago, a third of the world’s population became infected with the virus within the space of just a few weeks. The number of victims who succumbed to the Spanish flu far exceeds the death toll of the First World War. The book draws unmistakable parallels with the coronavirus pandemic as a global phenomenon and highlights the fragility of the civilised world even in the 21st century. The pillars on which our society is built – politics, the economy, social traditions, culture and education – have all been shaken to the core.

The fact that we aren’t completely helpless in the face of this virus – that we can defend ourselves and contain a pandemic caused by the SARS-CoV-2 coronavirus (severe acute respiratory syndrome coronavirus type 2) as a cause of COVID-19 – is largely thanks to the scientific progress we have made over the past few decades. Coronaviruses are widespread among mammals and birds around the world. They can expand their spectrum of hosts and cross over into different species with relative ease<sup>1</sup>. SARS-CoV-2 is an enveloped RNA virus that forms virions with large surface proteins (spike proteins) and has a single-stranded RNA genome.

Effective tests for all groups in the population are at the centre of the fight against SARS-CoV-2, and Germany has rolled out a comprehensive testing programme. Targeted testing allows scientists to rapidly and accurately determine the number of infections in Germany and their distribution across the country, helping to provide a more up-to-date and more accurate picture of the current situation. This information is essential to break chains of infection and to prevent our healthcare system from being overwhelmed.<sup>2</sup>

“We want to nip this virus in the bud. We can only achieve this with preventive mass testing in hospitals and care homes and by testing all of the contacts of an infected person where it is possible to do so. Money

should be no object. It is much more costly to not test enough than to test too much.” (Jens Spahn). Germany has developed a national testing strategy (see figure 1), but various organisations have also implemented their own testing concepts. Since 8 March 2021, the German federal government has covered the costs of regular rapid antigen tests for all citizens. Tests can be obtained from official state or local authority test centres and stations or from other authorised test points (such as doctors’ surgeries or pharmacies). To deliver this programme, the government has secured a share of up to 925 million rapid antigen tests. Rapid antigen tests must be carried out by trained staff.<sup>3</sup>

Depending on the clinical situation and question, if there is reason to suspect that a person is infected with the novel coronavirus SARS-CoV-2, investigative material should be taken from the upper respiratory tract. Where possible and if clinically required, samples should also be taken from the lower respiratory tract. To enter cells, SARS-CoV-2 uses receptors

**Author**  
 Dr Sabine Kaufmann  
 Head of CSSD  
 Klinikum Saarbrücken gGmbH  
 Winterberg 1  
 66119 Saarbrücken, Germany  
 www.klinikum-saarbruecken.de


Fig. 1: National testing strategy for SARS-CoV-2.

and enzymes that are co-expressed at a high level in the olfactory epithelium; scientists believe that this is how SARS-CoV-2 multiplies so effectively in, and is expelled so easily from, the upper respiratory tract.<sup>4</sup>

*Upper respiratory tract:*

- Nasopharynx swab (nose/throat swab)
- Oropharynx swab (throat swab)

*Lower respiratory tract:*

- Bronchoalveolar lavage
- Sputum (in patients with a productive cough; PPE required)
- Tracheal secretions

Nasopharynx swabs are the gold standard when taking samples from the upper respiratory tract for SARS-CoV-2 testing<sup>5</sup>. Throat swabbing is usually more tolerable than this type of swabbing for most patients; the molecular diagnostic sensitivity is comparable<sup>6</sup> or slightly reduced<sup>7</sup>. Throat and nose swabs can also be combined. The sampling stage is critical to the results of the test.

Since March 2020, testing capacity to check for the presence of genetic material of the virus (PCR testing) has been continually expanded; at the time of writing, one year later, over two million PCR tests can be performed each week. But what does PCR – a term we've all heard countless times since the start of the pandemic – actually mean?

PCR stands for polymerase chain reaction. The technique – invented by biochemist Kary Mullis, who won the 1993 Nobel Prize in Chemistry in recognition of his work – allows very small volumes of genetic material of a pathogen to be amplified or copied directly in vitro directly from a patient sample. The method only copies a short, precisely defined part of the genome, a single gene or part of a gene. This 'copying' stage is what makes PCR tests different from rapid antigen tests, which do not copy the genetic material prior to testing. In addition to the genetic material of the section to be copied (the template), a PCR test requires a number of other 'ingredients', including specific primers to determine the start point of synthesis on the DNA. These primers – which are short nucleotide sequences – are complementary to DNA; each corresponds to a very

small section that it can specifically attach itself to. The primers set the start and end points of the range to be copied. The DNA polymerase – an enzyme that combines the components of DNA then replicates or synthesises them in the desired range – then steps in. The components of DNA are adenine, guanine, thymine and cytosine. A buffer solution provides suitable chemical conditions for the reaction.

The PCR takes place in a number of stages. During the first stage of denaturation, the complementary double strands of the DNA are separated. The DNA is turned into single strands. During the second stage, which is annealing, the temperature is reduced and the primers are activated. After annealing, the third stage – extension of the DNA – takes place. Attaching to the primers, the polymerase builds a new strand on each of the free strands of the template DNA. New double strands form. The completion of these steps – denaturing, annealing and extension – marks the end of the first cycle of the PCR. Each additional cycle doubles the amount of template DNA. In the case of SARS-CoV-2, the virus genome is comprised not of DNA, but of RNA. The first step is to convert the RNA into DNA using the enzyme reverse transcriptase. Once this has been done, the cycles proceed in the same way.

In the gold-standard SARS-CoV-2 PCR test, the type of reaction carried out is a so-called 'real-time PCR'. A real-time PCR (qPCR) allows scientists to monitor the progress of the polymerase chain reaction as it happens, in real time. While the PCR is taking place, the tester can see whether the defined section of a gene is present in the patient sample. In addition to gene-specific primers, fluorescent tags are used. These tags are also gene-specific and will only bind to the targeted gene section. During the extension phase carried out by the polymerase, a measurable fluorescent signal is emitted. The fluorescence increases in proportion to the amount of PCR product and can be used in quantification. The higher the initial amount of target nucleic acids – or in other words, virus – the earlier a significant increase in fluorescence is observed. The higher the starting concentration of the target sequence, the lower the number of temperature cycles required for fluorescence to rise. The point at which the fluorescence becomes linear is known as the Ct value (Ct = crossing threshold). The

value represents the number of temperature cycles up to this point. The Ct value therefore shows a clear correlation with the virus concentration in the investigative material.

To ensure quality in molecular diagnostics – regardless of whether the test being performed is a PCR or rapid antigen test – it is essential to perform quality checks such as positive and negative controls on all tests. These checks allow scientists to demonstrate the reproducibility of the tests based on the generated measured values; they also enable relevant parameters, such as the detection limit, to be defined, and any deviations from the expected test performance to be identified. The Ct value from the real-time PCR is only a semi-quantitative value, and values obtained in different laboratories are not directly comparable if there is no reference value. A precisely quantified standard can be used to convert the Ct values to an RNA copy number per reaction and, if necessary, per sample volume (source: Robert Koch Institute guidance on testing patients for infection with the novel coronavirus SARS-CoV-2).

As a general rule of thumb, there is no one test that is completely free of errors. The PCR tests used to detect SARS-CoV-2 all over the world are not all (equally) reliable, even in defined laboratory conditions. Their performance depends on a number of factors, including the virus target and the amount of virus RNA in the sample.

Real-time PCR tests are highly sensitive with a high specificity, which means that they can be affected by how the samples were taken, stored and transported. However, the risk of false results is still low and is further minimised by various quality checks. The Robert Koch Institute requires two or three different gene regions to be examined. This means that the tests show evidence of different viral RNA sections, e.g. coding for coating proteins, the nucleocapsid or the RNA-dependent RNA polymerase (RdRp). A test is only deemed positive when positive results are found across the board.<sup>8</sup>

In October 2020, the national testing strategy was extended to include the use of tests that show evidence of virus proteins. Depending on their design, these antigen tests are suitable for use at the patient's loca-

tion (rapid antigen test, or point-of-care/PoC antigen test), as a single test or as a laboratory test for investigating larger sample volumes. Antigen tests involving sampling by trained staff help to further expand testing capacity.

A swab is used to take a sample from the nose-throat area. The swab is placed in a buffer solution to dissolve the sample material. The dissolved sample material is then applied to a test strip in a test cassette using a dropper, together with a buffer solution. The test strip reacts to the presence of certain viral protein components by changing colour, in a similar way to a pregnancy test. The test strip also has a control line, which confirms that the test has been performed correctly. The result appears after 15 to 30 minutes.

Since the end of February 2021, citizens have also been able to access antigen tests for at-home use (self-tests). These tests are deployed as part of testing strategies in organisations such as schools/childcare facilities and companies,<sup>9</sup> and are available from pharmacies, health and beauty retailers and supermarkets. All tests for home use must be designed and manufactured to ensure that the safety and performance of the medical device (incl. instructions, labelling etc.) are guaranteed when the test is used by persons who are not health-care professionals. Result quality must also be guaranteed when the test is used under these conditions. This requirement includes the use of the test itself and a user-friendly, reliable process for sampling and indication of the results. The correct use of the test has a significant impact on the reliability of the result. Studies have shown that, with correct instruction, the samples and the associated antigen test results obtained by individuals are comparable with the samples and results obtained by trained medical staff. High-quality self-tests have a clearly legible label on the packaging stating that the product has been granted special approval by the German Federal Institute for Drugs and Medical Devices and/or a CE mark, together with four-digit reference number for the notified body. (Source: Robert Koch Institute, National testing strategy – who is tested for SARS-CoV-2 infection in Germany?).

Due to the low sensitivity and specificity of antigen tests, they are only a useful addition to other measures

under certain conditions. For an antigen test to show a positive result, a larger amount of virus must be present than for a PCR test, as the target is not multiplied before detection takes place (lower sensitivity). Antigen tests therefore primarily indicate a high viral load, which means that a negative antigen test result does not rule out the possibility of a SARS-CoV-2 infection. For this reason, these tests should only be used in situations where a false negative result would not have serious consequences (such as a failure to recognise an infection on admission to hospital). The result is heavily dependent on the time the sample was taken, the quality of the sample (e.g. nose swab) and whether or not the test has been performed correctly. Particularly in cases where the time of infection is not known (such as asymptomatic persons) and in the first seven days after infection, the viral load in the upper respiratory tract changes very quickly. A negative result on day 4 after infection could turn into a positive on a second test just one day later, as the virus will have had a chance to multiply in the nasopharynx. Furthermore, rapid antigen tests are not as highly specific as PCR tests. Rapid antigen tests can show a positive result even if the person is not infected, which is not the case with PCR tests performed in a laboratory. For this reason, a positive result on an antigen test should generally be confirmed with a PCR test.

Like rapid tests, self-tests provide an indication of what is going on at a specific moment in time. For practical reasons, the result can only be considered “valid” for a maximum of 24 hours. If testing is carried out repeatedly, there is a higher chance of picking up a transmissible infection at an early stage.

The German Federal Institute for Drugs and Medical Devices (BfArM) and the Paul Ehrlich Institute provide information on SARS-CoV-2 testing systems. The information provided by the two organisations is complementary; the legal framework is provided by a number of regulations, including the German Coronavirus Testing Ordinance (TestV). The BfArM provides a list of antigen tests that can be used to directly test for the presence of coronavirus SARS-CoV-2 pathogens.

The list includes tests which, according to information provided by the manufacturer, are legally marketed in Germany and Europe in accordance with the Medical Device Directive and meet all of the minimum criteria for antigen tests as defined by the Paul Ehrlich Institute in consultation with the Robert Koch Institute (RKI). The list is updated on an ongoing basis. (Source: [www.pei.de](http://www.pei.de)).

#### References

1. Graham and Baric, 2010.
2. Source: Robert Koch Institute, National testing strategy – who is tested for SARS-CoV-2 infection in Germany?
3. [www.bundesgesundheitsministerium.de](http://www.bundesgesundheitsministerium.de).
4. Sungnak et al., 2020). (Source: Robert Koch Institute, SARS-CoV-2: Virological basic data and virus variants, as at 07.04.2021.
5. WHO, 2020b.
6. Wolfel et al., 2020.
7. Covid-Investigation Team, 2020; Wang et al., 2020.
8. Source: [aerzteblatt.de](http://aerzteblatt.de), Medizinreport, PCR tests for SARS-CoV-2: Interpreting results correctly.
9. Source: Robert Koch Institute, National testing strategy – who is tested for SARS-CoV-2 infection in Germany?


## New Instrument Reprocessing Working Group brochure in the series: “Instrument Reprocessing: Reprocessing of Instruments to Retain Value”

Forty-five years ago, a group of experts founded the Instrument Reprocessing Working Group (AKI), with the aim of gathering, bundling and publishing knowledge from the specialist areas of medical device development and production, washer-disinfectors and sterilisers, as well as process chemicals and their interactions during reprocessing. Based on scientific findings and experience in the reprocessing of reusable medical devices, the group compiles and publishes practical guidance with a focus on preventive value retention. Since its release in 1979, almost 400,000 copies of the group’s first publication – the “Red Brochure” entitled “Instrument Reprocessing: Reprocessing of Instruments to Retain Value” – have been circulated around the world, with its advice translated into 20 languages. The “Red Brochure” is an important point of reference for users in many countries, as well as for employees in central sterile services departments (CSSDs).

One of the key chapters in the “Red Brochure” tackles surface changes and damage to instruments observed during reprocessing. The guidance sets out the possible causes and risks, provides recommendations as to how to avoid this kind of damage and explains whether changes can be reversed using appropriate measures or whether repair is required.

During a 2019 European conference on endoscope reprocessing, one of the sessions looked at the damage that typically occurs to an endoscope during reprocessing. It was found that, in many cases, the employees in the RUMED (Reprocessing Unit for Medical Devices) did not have sufficient knowledge of how to identify, interpret and assess damage to flexible endoscopes. During the discussion, participants came up with the idea of creating a document similar to the “Red Brochure” to provide guidance on reprocessing flexible endoscopes in a way that would retain value.

**Fig. 1:** New brochure for Instrument Reprocessing.

With expert support from leading endoscope manufacturers, the AKI has now made this idea a reality. The new brochure in the “Reprocessing of Instruments to Retain Value” series focuses specifically on value retention in the reprocessing of flexible endoscopes. The brochure aims to provide reprocessing personnel with recommendations and guidance for the correct and safe reprocessing of flexible endoscopes. In turn, this will help to retain the function and value of the devices for long periods of time.

The first edition of the brochure is set to be unveiled at the WFHSS Congress in November 2021 and will be available for download from [www.a-k-i.org](http://www.a-k-i.org) after publication.

# Efficacy of airborne room disinfection – European test standard EN 17272<sup>1</sup>

## Authors

Henrik Gabriel  
Laboratory manager  
Dr. Brill + Partner GmbH  
Institute for Hygiene and Microbiology  
Stiegstück 34, 22339 Hamburg, Germany  
henrik.gabriel@brillhygiene.com  
www.brillhygiene.com

Dr Dajana Paulmann  
Deputy Laboratory Manager  
Dr. Brill + Partner GmbH  
Institute for Hygiene and Microbiology  
Norderoog 2, 28259 Bremen, Germany  
dajana.paulmann@brillhygiene.com  
www.brillhygiene.com

Dr Florian H. H. Brill  
Executive Director and Co-Proprietor  
Dr. Brill + Partner GmbH  
Institute for Hygiene and Microbiology  
Norderoog 2, 28259 Bremen, Germany  
florian.b@brillhygiene.com  
www.brillhygiene.com

Henrik Gabriel, Dajana Paulmann,  
Florian H. H. Brill

For over 100 years, airborne processes have been deployed as an alternative and complement to manual surface disinfection methods involving spraying or wiping<sup>2</sup>. The advantage of these airborne methods is that they can potentially reach all of the surfaces in a room. Although formaldehyde vaporisation was the dominant method in Germany for many years, more recent developments have seen new active substances and techniques added to the mix, including methods that use hydrogen peroxide, peracetic acid or ozone. Just like conventional disinfection procedures, all of these methods need to be tested for efficacy. However, the established methods for testing chemical surface disinfectants did not adequately reflect how these products were being used in practice, and there were no harmonised, practical European standards.

Generally, these kinds of processes are used as complements to “standard” surface cleaning and disinfection, and they cannot replace these methods. They may be used in areas such as infection rooms and operating theatres, where airborne disinfection methods can be a useful addition to the disinfection process.

When chemical disinfectants are atomised, vaporised or dispensed via an aerosol, the solutions used can be tested in suspension tests such as those outlined in EN 13727<sup>3</sup> and EN 14476<sup>4</sup>, or using existing practical methods such as the four-field test (EN 16615)<sup>5</sup>. In principle, gaseous agents, such as ozone, cannot be tested using this procedure. Even if testing is technically feasible, the application conditions for airborne room disinfection are so unique that the test is not expected to yield valid results. The test may underestimate or overestimate the efficacy of the procedure.

In the early 1980s, in direct response to this issue, a practical testing method was established in France in the form of NF T 72-281<sup>6</sup>. This method has been refined over the past few decades and was published as DIN EN 17272 in June 2020.

Intended for use in different room sizes

The method uses “small” and “large” rooms to test two common areas of application. Small rooms reflect use in incubators or cold storage cells, while large rooms represent production rooms or operating theatres. In the laboratory tests, rooms between 0.25 m<sup>3</sup> and 4.00 m<sup>3</sup> represent small rooms, while rooms between 30 m<sup>3</sup> and 150 m<sup>3</sup> are used for large rooms.

As the active substances react on surfaces or at least coat them, the test rooms should not be furnished. They must be capable of being sealed to prevent the entry of gases to prevent the loss of the active substance through diffusion. They must also maintain a consistent temperature (20 °C ± 2 °C) and relative humidity (50 % to 75 %). The rooms must be well ventilated so that any remaining traces of the active substance can be allowed to escape once testing is complete.

In practice, however, test rooms do not satisfy these test room requirements, so test results obtained in accordance with DIN EN 17272 are not directly transferable to real-life applications. The aim of the testing process is to obtain reproducible results to assess the general efficacy of the test and to allow for comparison between different application conditions and procedures. For this reason, the standard recommends that additional checks are performed to verify the suitability of the test for the specific local application conditions.

## Testing principle

To determine efficacy, the test follows a practical model in which an application device is used to distribute the active substances in the room. The test looks at either a combination of a device and an active substance solution, or a generator that generates the active substance directly.

A mixture of test organisms/viruses and an organic load is applied to a stainless steel test specimen, dried and then exposed to the disinfection process. The number of surviving test organisms/infectious virus particles is then determined and compared with the untreated control variant. The resulting logarithmic reduction factor (Rf) reflects the decline in living test organisms and therefore the efficacy of the treatment. In the case of viruses, the decline in infectiousness is measured using the TCID50 process.

The standard specifies mandatory test conditions, but these conditions may also be modified by manufacturer instructions. Variables such as the selected test organisms, the process temperature and humidity levels are all subject to change. The aim is to ensure that the test conditions reflect real-life use as far as possible. In terms of test organism selection, the standard builds on the proven substitution system to illustrate the effect on various groups of organisms. Unlike other standards, it sets out combined test guidelines for various test organisms and viruses in a single standard. DIN EN 17272 is the only standard to cover the assessment of bactericidal, yeasticidal/fungicidal, sporicidal, mycobactericidal and virucidal activity (against human and animal pathogenic viruses and bacteriophages) from the test areas of human medicine, veterinary medicine and food, industry, domestic and public facilities. Table 1 lists the mandatory test organisms and the required reductions.

The organic load substances used are known combinations described in other standards and specific to the test areas. These combinations are listed in table 2. For test organisms where drying is critical – such as *Pseudomonas aeruginosa*, *Escherichia coli*, *Proteus hauseri* or *Candida albicans* – 0.5 % skimmed milk can be added to each load to prevent drying. Stainless steel plates are used as the test surface. Measuring

3–4 cm in diameter, these plates are slightly larger than those used in DIN EN 13697<sup>7</sup> and EN 16777<sup>8</sup>, which have a diameter of 2 cm. As a result, 50 µl of the test organism and load mix can be distributed across an area of around 3 cm<sup>2</sup>, which is roughly equal to the entire surface of the fomite used up to this point. This produces a thinner layer, which in turn means that the test organisms that survive drying are exposed to the influence of the test procedure equally. In EN 13697, for example, the suspension is applied in drops, which means that the layer thickness in the centre is significantly greater than at the edge. Depending on the room volume, three process challenge devices for each process and test organism are affixed at defined heights and distances so that the inoculated side is facing away from the release source. Figure 1 shows a diagram of a test room setup. Two simultaneously created control process challenge devices are not exposed to the test, but are stored in the laboratory for the duration of the reaction time in temperature and relative humidity conditions that are as close to the test conditions as possible. In procedures with high humidity and long reaction times (usually over an hour), this means that bacteria sensitive to drying are significantly more likely to survive. As is the case in other standards designed to test efficacy, the validity of the test is ensured using controls. As well as determining compliance with the prescribed bacterial count range

Effective range	Test organisms	Test area					
		Human		Veterinary		Industry	
		Selected	Rf	Selected	Rf	Selected	Rf
Bactericidal	<i>Staphylococcus aureus</i>	X	5	X	5	X	5
	<i>Enterococcus hirae</i>	X	5	X	5	X	5
	<i>Escherichia coli</i>	X	5	-	-	X	5
	<i>Pseudomonas aeruginosa</i>	-	-	X	5	X	5
	<i>Acinetobacter baumannii</i>	X	5	-	-	-	-
	<i>Proteus hauseri</i>	-	-	X	5	-	-
Fungicidal	<i>Candida albicans</i>	X	4	X	4	X	4
	<i>Aspergillus brasiliensis</i>	X	4	X	4	X	4
Yeasticidal	<i>Candida albicans</i>	X	4	X	4	X	4
Sporicidal	<i>Bacillus subtilis</i>	X	4	X	3	X	3
Mycobactericidal	<i>Mycobacterium terrae</i>	X	4	-	-	X	4
	<i>Mycobacterium avium</i>	X	4	X	4	X	4
Tuberculocidal	<i>Mycobacterium terrae</i>	X	4	-	-	X	4
Virucidal	Murine norovirus	X	4	-	-	X	4
	Adenovirus, type 5	X	4	-	-	X	4
	Porcine parvovirus	-	-	X	4	-	-
Phage effect	<i>Lactobacillus lactis</i> P001	-	-	-	-	X	4
	<i>Lactobacillus lactis</i> P008	-	-	-	-	X	4

**Tab. 1:** Selection of test organisms and required reductions for various test areas (Selected = organism selection).

Test area		
Human	Veterinary	Industry
Mandatory conditions		
Low load: 0.3 g/L BSA	High load (low level): 3.0 g/L BSA	Low load: 0.3 g/L BSA
Additional conditions		
High load: 3.0 g/L BSA + 3.0 ml/L sheep erythrocytes	High load (high level): 10.0 g/L BSA + 10 g/L yeast extract	High load: 3.0 g/L BSA

**Tab. 2:** Selection of organic loads for test areas (BSA = bovine serum albumin).

and the previously mentioned water control, the test also checks that the inactivation procedure is sufficiently effective. To do so, process challenge devices (with the same organic load as the test) are exposed to the process, and the toxicity of the recovery solution and the process challenge device are determined after recovery.

### Reaction time?

In practice, a distinction is made between the time in which the effective concentration is built up and the time after this point until the end of the test as the reaction time. In laboratory tests, the time at which the process challenge device is removed marks the end of the reaction time. With appropriate PPE, the device can be removed as soon as the process time specified by the manufacturer has elapsed or after a decontamination phase specified by the device.

In the standard, the reaction time is the entire process time from the initial release of the product up to the point at which the process challenge device is removed. This is a logical approach, because the two phases cannot be separated in terms of their effect on reduction. A disinfection procedure with a short start phase and the same overall duration will not necessarily display similar levels of efficacy, even if the same amount of active substance were to be distributed within a shorter period of time. This must be checked separately.

Under the mandatory conditions of the standard, the reaction time must be less than 15 hours. However, this is generally a challenge for many test organisms.

### New challenge – the distribution test

A new aspect of EN 17272 is the distribution test. In this test, further process challenge devices inoculated with *Staphylococcus aureus* are placed in the corners of the test room. Two are placed in opposite corners of the room, spaced 15 cm (in small rooms) and 50 cm (in large rooms) from the floor; another two are secured on the ceiling. In each corner, one

process challenge device is positioned vertically, facing away from the source, while another is placed horizontally, facing the ceiling (if mounted on the ceiling) or the floor.

In test conditions, the tested procedure must achieve a reduction of at least five log levels. The distribution test can be performed at the same time as the actual test or as a preliminary test under identical test conditions.

In practice, this test has been somewhat of a stumbling block. Many procedures that offer sufficient reduction for standard process challenge devices fail the distribution test, falling one or more levels short of the required five log levels. For this reason, the test can be a good preliminary test in many cases.

### Outlook

During the coronavirus pandemic, interest in automated room disinfection procedures has increased significantly. In an automated room disinfection process with ozone as the active substance, for example, tests demonstrated efficacy against the bacteriophage  $\Phi 6$  ( $\phi 6$ )<sup>10</sup> and the bovine coronavirus as a surrogate for the SARS-CoV-2 virus<sup>9,10</sup>. This procedure uses chemical processes as well as processes that rely on physical principles (e.g. UV-C radiation). It is important to bear in mind that although EN 17272 describes a practical test procedure, it addresses a process based on a chemical effect. The standard test organisms were selected based on their known high resistance to chemical disinfectants. For other methods, such as UV radiation, other practice-relevant organisms may be more suitable for testing. Shading also plays a larger role than with traditional organic loads. In these applications, modifications of the method should be considered.

On the other hand, for airborne chemical disinfection procedures, EN 17272 is a highly practical method, carried out in a similar way to a controlled field test. As known process challenge devices and standard organisms are used, test laboratories experienced in testing surface disinfectants find this method easy to implement.

However, these laboratories do need to have access to suitable test rooms, which can be a major hurdle in many cases. Not all facilities have access to a room that is at least 30 m<sup>3</sup> with sufficient sealing capacity. However, for testing in smaller rooms, existing equipment (e.g. refrigerators) can be repurposed or suitable test boxes can be constructed with relatively little effort.

Initial experiences suggest that the new distribution test may be a major obstacle for many conventional processes. The problem seems to be that the active substance does not travel evenly into the corners of the room. Furthermore, the long reaction times mean that the issues with drying-sensitive test organisms observed with other methods become even more prevalent. If 0.5 % skimmed milk were to be added to the standard loads for all drying-sensitive test organisms for practical reasons, the additional organic load alone would result in a significant decline in efficacy.

It is also important to note that the test room represents an ideal set of conditions. In practice, it is very difficult to find an unfurnished room with completely inert surfaces. With this in mind, the results of DIN EN 17272 cannot be directly transferred to practice; tests in sample rooms at the location of use are recommended to ensure that the substance will be effective for the user.

In spite of these limitations, DIN EN 17272 is a good way to standardise efficacy testing for airborne disinfection procedures and increases the safety of disinfection as an infection prevention measure.


Fig. 1: Diagram of a test room.

## References

1. German version, EN 17272. Chemical disinfectants and antiseptics. Methods of airborne room disinfection by automated process. Determination of bactericidal, mycobactericidal, sporicidal, fungicidal, yeasticidal, virucidal and phagocidal activities. DIN EN 17272:2020-06.
2. Koch, R. Über Desinfektion. Mitteilungen aus dem Kaiserl. Gesundheitsamte, 1881, Bd. I, Berlin.
3. German version, EN 13727. Chemical disinfectants and antiseptics. Quantitative suspension test for the evaluation of bactericidal activity in the medical area. Test method and requirements (phase 2, step 1). DIN EN 13727:2012+A2:2015.
4. German version, EN 14476. Chemical disinfectants and antiseptics. Quantitative suspension test for the evaluation of virucidal activity in the medical area. Test method and requirements (phase 2/step 1). DIN EN 14476:2013+A2:2019.
5. German version, EN 16615. Chemical disinfectants and antiseptics. Quantitative test method for the evaluation of bactericidal and yeasticidal activity on non-porous surfaces with mechanical action employing wipes in the medical area. Test method and requirements (phase 2, step 2). DIN EN 16615:2015-06.
6. Norme française NF T 72-281. Methods of airborne disinfection of surfaces – Determination of bactericidal, fungicidal, yeasticidal, mycobactericidal, tuberculocidal, sporicidal and virucidal activity, including bacteriophages. NF T 72-281 8 November 2014.
7. German version, EN 13697. Chemical disinfectants and antiseptics, Quantitative non-porous surface test for the evaluation of bactericidal and/or fungicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas. Test method and requirements without mechanical action (phase 2, step 2). DIN EN 13697:2015+A1:2019.
8. German version, EN 16777. Chemical disinfectants and antiseptics. Quantitative non-porous surface test without mechanical action for the evaluation of virucidal activity of chemical disinfectants used in the medical area. Test method and requirements (phase 2/step 2). DIN EN 16777:2018.
9. Brill, F H H et al. Virucidal efficacy of an ozone-generating system for automated room disinfection. J Hosp Infect. 2021 June 14, doi: 10.1016/j.jhin.2021.06.004.
10. Knobloch, J K. An automated room disinfection system using ozone is highly active against surrogates for SARS-CoV-2. J Hosp Infect. 2021 Jun; 112:108-113, doi: 10.1016/j.jhin.2021.04.007.

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## Introduction

During the last millennium, Clostridioides difficile infection (CDI, previously known as Clostridium difficile) was a rare and generally harmless disease with a low specific mortality rate,<sup>1</sup> observed only in severely ill patients following antibiotic treatment and in newborn and premature babies.

However, after the year 2000, rates of CDI began to climb, initially in North America.<sup>1</sup> Cases were seen in previously healthy people, and infections were proving serious and even fatal. Early examples included a 50-year-old male who had undergone elective hip replacement surgery<sup>1</sup> and fatal outcomes for pregnant women.<sup>2</sup> The new epidemic strains responsible for these developments were found to possess additional virulence factors, which were driving the significantly higher rates of disease and mortality.<sup>1</sup>

In 2003, these new strains reached Europe.

The first outbreaks were identified in Stoke Mandeville Hospital in the UK, which recorded a total of 498 CDI cases between 2003 and 2006, 127 of which were fatal.<sup>3</sup>

Germany's first cases – which were also directly associated with an outbreak – were identified in Trier in 2007. Subsequent investigations by local authorities and the Robert Koch Institute (RKI) showed that the disease was already endemic in the region by this point.<sup>4</sup> Today, Clostridioides difficile is one of the most common culprits behind nosocomial infections. The CDC even lists Clostridioides difficile as one of the five most urgent threats to health.<sup>5</sup> This decision was based on data from 2017, a year in which almost 13,000 people died from Clostridioides difficile infections in the USA and which saw a quarter of a million cases in US hospitals, resulting in costs of 1 billion dollars for the healthcare sector.<sup>5</sup>

In the European Point Prevalence Survey (PPS) carried out by the ECDC in 2016/2017, CDI was found to ac-

count for almost 5 % of all nosocomial infections in the EU.<sup>6</sup> Data gathered from Germany in 2016 shows that CDI was behind 10 % of these infections.<sup>7</sup>

In 2016, Cassini et al. calculated the impact of the six main nosocomial infections in the EU. CDI was in fifth place with over 150,000 cases per year and 31.2 disability-adjusted life years (DALYs) per 100.00 citizens (equivalent to 6 % of DALYs overall).<sup>8</sup>

## CDI surveillance at mkm

Marienhaus Klinikum Mainz (mkm) is categorised as a Level II hospital under the German system and has 602 beds and around 1,500 employees. Each year, the hospital provides in-patient and outpatient care for 50,000 people across 19 clinics and 10 centres.<sup>9</sup> Since 1999, mkm has been carrying out surveillance of nosocomial infections via KISS, a national infection surveillance system operated by the German National Reference Center for Surveillance of Nosocomial Infections (NRZ).

When Hubert Holz heard and read media reports about the Stoke Mandeville outbreak while on holiday in the UK in 2005, he decided to establish his own surveillance system for CDI at mkm, with the aim of identifying any upward trend in case numbers at an early stage so that action could be taken as necessary.

It soon became clear what a good idea this was. Once the plan had been approved by the hygiene committee and the medical director, mkm embarked on its own surveillance programme for Clostridioides difficile on 1 January 2006. In the first year of the programme, 100 patients with clinical manifestations of and microbiologically confirmed CDI were recorded (see fig. 1). No-one had expected to see such high case numbers. However, the first new epidemic strain (then known as RT 027) was only detected at mkm at the end of 2007.

From 2007 onwards, the NRZ also offered its own CDAD-KISS module for recording CDI; initially, 34 other German hospitals participated in this scheme alongside mkm<sup>10</sup> (but by 2018, this number had risen to 571!).<sup>11</sup>

### Number of CDI-cases mkm 2006 - 2020


Fig. 1: Development of CDI cases at mkm, 2006–2020.

### Development of case numbers at mkm and in Germany

Over the next few years, case numbers rose significantly, not only at mkm but also in Germany (see fig. 2) and across the world. As well as embarking on a targeted prospective surveillance programme, case numbers prompted mkm to implement two further measures to combat CDI as early as 2006. The hospital developed an information sheet setting out hygiene measures specific to CDI and provided laminated and electronic copies for all wards.<sup>12</sup> 2006 also saw the development of the first antibiotic therapy

guidelines for mkm, which laid the foundations for the Antibiotic Stewardship system in the organisation.<sup>13</sup>

In 2008, mkm also introduced the hygiene consultation concept for patients with CDI.<sup>14</sup> In this system, an on-site hygiene specialist discusses the patient’s stay and care at the ward with the responsible members of staff. This means that the required additional hygiene and barrier measures (such as sporicidal surface disinfection and additional hand-washing) can be immediately communicated to the staff responsible for the patient’s care and any gaps in their knowledge can be filled. From 2009 to 2015,

### Comparison of overall prevalence (CDI-cases / 100 patients) mkm and standardized infection ratio NRZ


Fig. 2: Comparison of overall prevalence of CDI cases, mkm and NRZ 11 2006–2020.

CDI case numbers remained stable at a very high level (see fig. 2).<sup>11</sup> During this period, the government also set out a number of measures in an attempt to contain CDI, including an obligation to monitor the use of antibiotics (from 2011), an obligation to conduct CDI surveillance in all hospitals (from 2013) and the initiation of the Antibiotic Consumption Surveillance (AVS) module by the RKI. These measures started to pay off in 2016: Since then, both the prevalence of CDI on admission to hospital and the incidence of cases acquired in Germany<sup>11</sup> has been gradually declining (see fig. 2).

## Conclusion

15 years of surveillance of *Clostridioides difficile* infections at mkm prove how important it is to identify problematic nosocomial pathogens at an early stage. It is only possible to implement a tailored, targeted strategy to contain and combat infection, and to evaluate the efficacy of this strategy, if you have access to valid data. Alongside measures to isolate infected patients, using sporicidal disinfectants on surfaces and equipment and

stepping up hand hygiene practices, the strategy also involves the careful use of antibiotics and the introduction of a full Antibiotic Stewardship system. But we cannot eliminate the use of antibiotics completely! We still need to use antibiotics to fight infection, even if doing so could lead to CDIs. It is essential that we weigh up the risks and benefits of antibiotic therapy in each individual case. As antibiotics will remain a crucial part of our infection-fighting strategy, we do not believe that we will be able to eradicate *Clostridioides difficile* infections completely.

However, thanks to the measures that have been implemented in mkm and across Germany as a whole, CDI case numbers are now almost down to the levels we saw many years ago, before the epidemic strain began to spread around the world. This is clear evidence of how preventive measures can be used to contain and manage major epidemic events; something that is a great source of hope in the current SARS-CoV-2 pandemic.

## References

1. Clostridium-difficile-Infektionen: Nosokomiale Ausbrüche durch einen neuen, besonders virulenten Stamm in den USA, Kanada, England, Belgien, Holland und Frankreich; Kist, Manfred et al.; Epidemiologisches Bulletin 2006; 36 (11); 309-311.
2. Clostridium difficile-associated diarrhea: an emerging threat to pregnant women; Roupheal, Nadine G. et al.; Am J Obstet Gynecol 2008 (98) :635.e1-635.e6.; <https://doi.org/10.1016/j.ajog.2008.01.062>.
3. Investigation into outbreaks of Clostridium difficile at Stoke Mandeville Hospital, Buckinghamshire Hospitals NHS Trust; Commission for Healthcare Audit and Inspection; 2006; ISBN 1-84562-103-4.
4. Clostridium-difficile-Ribotyp 027: Epidemiologie und Klinik des erstmaligen endemischen Auftretens in Deutschland; Jansen, Andreas et al.; Z Gastroenterol 2010; 48(9): 1120-1125 ; <https://doi.org/10.1055/s-0029-1245269>.
5. Antibiotic Resistance Threats in the United States; U.S. Department of Health and Human Services; 2019; [www.cdc.gov/DrugResistance/Biggest-Threats.html](http://www.cdc.gov/DrugResistance/Biggest-Threats.html). <http://dx.doi.org/10.15620/cdc:82532>.
6. Prevalence of healthcare-associated infections, estimated incidence and composite antimicrobial resistance index in acute care hospitals and long-term care facilities: results from two European point prevalence surveys, 2016 to 2017; Suetens, Carl et al.; Euro Surveill. 2018; 23 (46): 1-18; <https://doi.org/10.2807/1560-7917.ES.2018.23.46.1800516>.
7. Prävalenz von nosokomialen Infektionen und Antibiotika-Anwendung in deutschen Krankenhäusern; Behnke, Michael et al.; Dtsch Arztebl Int 2017; (114): 851-7. <https://doi.org/10.3238/arztebl.2017.0851>.
8. Attributable deaths and disability-adjusted life years caused by infections with antibiotic-resistant bacteria in the EU and the European Economic Area in 2015: a population-level modelling analysis; Cassini Alessandro Cassini et al.; Lancet Infect Dis. 2019 Jan;19(1):56-66; <http://dx.doi.org/10.1016/j.laninf.2018.10.021>.
9. Menschlich und kompetent. [https://www.marienhaus-klinikum-mainz.de/fileadmin/user\\_upload/kkm\\_Imagebroschuere\\_2019.pdf](https://www.marienhaus-klinikum-mainz.de/fileadmin/user_upload/kkm_Imagebroschuere_2019.pdf).
10. Inzidenz der Clostridium difficile-assoziierten Diarrhö: Erste Ergebnisse aus dem CDAD-Modul des Krankenhaus-Infektions-Surveillance-Systems; Weitzel-Kage, Doris et al.; HygMed 2008; 33 (9): 353-356. <https://doi.org/10.25646/919>.
11. CDAD-KISS module – reference data; German National Reference Center for Surveillance of Nosocomial Infections; [https://www.nrz-hygiene.de/fileadmin/nrz/module/cdad/201901\\_201912-CDAD\\_Ref\\_Dev2.pdf](https://www.nrz-hygiene.de/fileadmin/nrz/module/cdad/201901_201912-CDAD_Ref_Dev2.pdf), version: 29 April 2019.
12. Hygiene management in in-patient care; Weiser, Tanja; Holz, Hubert; Aseptica 2014; 20 (1): 9-15.
13. Establishment of an Antibiotic Stewardship system based on the example of the Katholisches Klinikum in Mainz; Holz, Hubert; Aseptica 2019; 25 (1): 7-11.
14. The hygiene consultation concept at mkm; Holz, Hubert; Kiesel, Markus; Kiesel, Heike; Aseptica 2020; 26 (1): 9-13.

## New pressure/temperature data logger EBI 12-TP237

The new pressure/temperature data logger EBI 12-TP237 is suitable both for routine checks in CSSDs and for the validation of washer-disinfector, endoscopy washer-disinfector and steam sterilisation processes.


The data logger boasts a measurement range of 0–140 °C/1–4000 mbar with an accuracy of +/- 0.1 °C and +/- 20 mbar, which is confirmed on the calibration certificate supplied with the device. It comes with a practical hose connection which can be used to connect the device to the washer-disinfector/endoscopy washer-disinfector to measure the rinse pressure. The device can also measure the A0 value of the washer-disinfector at the same time as the rinse pressure. In steam sterilisers, the data logger records the temperature and pressure, enabling the theoretical saturated steam temperature to be calculated. The data logger can be used with the existing Xylem/ebro software systems Winlog.med and Winlog.validation once the free update to version 3.73 has been installed.

The data logger is equipped with sensors that comply with EN 285, EN 13060 and ISO 15883 and represents exceptional value for money.

Fig. 1: EBI 12-TP237 temperature/pressure data logger in use in a washer-disinfector.


Scan the code to visit the ebro shop:

## “Miele MOVE” portal monitors laundry machines

Has my programme completed properly and has disinfection taken place? When is it due to finish? The new “Miele MOVE” portal is designed to answer exactly these kinds of questions. The platform for PCs and mobile end devices simplifies a multitude of processes and provides a secure documentation solution. The platform is now available in Germany for the new Benchmark Machines, Miele’s latest generation of laundry machines, and will be expanded to many other appliances in the future. The Little Giants washing machines and tumble dryers are next in line, followed by a number of fresh water dishwashers. Miele MOVE will also be available in Austria and Switzerland from 2022.

# Requirements for the validation of cleaning and disinfection processes – DIN 58341

## Author

Iven Kruse  
Xylem Analytics Germany Sales  
GmbH & Co. KG, ebro  
Peringerstraße 10  
85055 Ingolstadt, Germany  
T: +49 841 95478-0  
F: +49 841 95478-80  
Iven.kruse@xylem.com

### Iven Kruse

The new DIN 58341 standard was published a year ago. Now is the perfect time to look back and ask: What has improved as a result of the new standard and what opportunities has it created?

DIN 58341 describes the requirements for the validation of cleaning and disinfection procedures. The standard provides a more precise description of the requirements of the existing DIN EN ISO 15883 standards, as well as the requirements of the guidelines published by the German Society for Hospital Hygiene (DGKH), the German Society for Sterile Supplies (DGSV) and the Instrument Reprocessing Working Group (AKI) for the validation and routine monitoring of machine-based cleaning and thermal disinfection processes; the guidelines published by DGKH, DGSV, the German Society for Digestive and Metabolic Disease (GDVS), the German Society of Endoscopy Nurses and Associates (DEGEA) and AKI for the validation and routine maintenance of machine-based cleaning and thermolabile endoscopy; and the guidelines published by DGKH, DGSV and AKI in partnership with the German Association for Applied Hygiene (VAH) on the validation of manual cleaning and manual chemical disinfection of medical devices.<sup>1,2,3,4,5</sup>

Under §8 of the German Medical Device Ordinance (MP-BetriebV), medical devices may only be used if they have been reprocessed using an appropriate validated procedure. The operator is responsible for ensuring that this is the case. The validation of the reprocessing procedure must be conducted on the instruction of the operator by qualified specialists/validators certified in accordance with §5. Compliance with the requirements set out in §5, para. 1 is demonstrated by the issue of a certificate by the responsible authority or another recognised body, such as ZLG, TÜV or a responsible authority in another EU country.

Certification must be implemented by 1 January 2020. The standard DIN 58341 is a tool to help meet the validation re-

quirements; it can be used to plan and conduct validations with the help of the applicable standards and guidelines.<sup>6</sup>

Section 5 of DIN 58341 describes the requirements for the performance of validations and the knowledge required to perform validations.

Under General, the standard describes the tests and assessments for the various process-relevant procedures in washer-disinfectors, in endoscopy washer-disinfectors and for manual cleaning and disinfection.

Furthermore, the validator must have knowledge of the medical devices to be reprocessed, the standards and guidelines, equipment and tools, testing and measuring devices, microbiology, chemistry, hygiene, devices and the process itself. For medical devices, knowledge of the design, risk assessment and product families is important. In the standards and guidelines, section 5.3 provides a more detailed description of the required knowledge of regulations, standards, directives, guidelines and recommendations. It is essential that the validator knows the current regulations and performs the validation in line with these regulations.

The validator must possess knowledge of the equipment to be used, including work tables, dosing technology, load carriers, adapters, pre-rinse devices, endoscopy fit testers, drying aids, compressed air guns, brushes, tim-

**Fig. 1: Washer-disinfector process validation.**

ers, water guns and water filters. The equipment and the DIN EN ISO 15883-compliant data logger are an important part of the validation process. The validator must also possess knowledge of the measurement technology, its accuracy, the design of the temperature sensor, the pressure sensor and its connections (e.g. Luer Lock) and the validation software used.

The validator must understand how the measurement technology is used during validation, how often it needs to be calibrated and how measurement uncertainty impacts on the validation result. To ensure that the complex validation software is operated correctly, the user must be trained in its functions. The standard requires microbiological knowledge to determine the quality of the process and rinse water; this knowledge is also important when using microbiologically contaminated process challenge devices and when determining the microbiological status of the reprocessed medical device. The test laboratory performing the tests must demonstrate knowledge of verifying and identifying microorganisms, pathogenicity, growth media, the temperature and duration of incubation and the determination and evaluation of colony numbers.<sup>1</sup>

In terms of cleaning tests, DIN 58341 requires knowledge of protein residue testing using the ortho-phthalaldehyde (OPA) method or the bicinchoninic acid (BCA) method, as well as photometry and potential interference factors. If necessary, negative controls must be carried out, regardless of whether the cleaning test checks are being performed on site or in a test laboratory. This ensures that the principle behind the method is sound and that interference factors that could lead to false-positive results are minimised and eliminated where possible.<sup>1,8</sup> During validation, the staff must comply with the operator hygiene standards set out in the hygiene plan. This includes wearing the appropriate clothing and personal protective equipment as well as wearing appropriate gloves during the sampling process.

Section 5.9 of DIN 58341 provides explanatory examples of the required machine and process knowledge. To perform a professional validation, the validator must possess electrotechnical, machine-specific, process-related and washer-disinfector/endoscopy washer-disinfector-specific knowledge. Electrotechnical knowledge encompasses the machine-specific safety equipment, the electrical functions such as cable breaks or short circuits and the electrotechnical layout of the machine, including the circuit diagram.

**Fig. 2: Pressure measurement during washer-disinfector process validation.**

For qualification, the validator must be familiar with safety functions such as the door locking mechanism, error messages, test functions and interrupting the program. Knowing how to interrupt the program is important as the program needs to be stopped at the end of the cleaning process to remove the cleaned process challenge device. Other important knowledge for the validator includes:

- Applications of the machine
- Load carriers and adapters
- The minimum and maximum rinse pressure and how to measure this pressure
- Defined rinse arm speed
- Type of heating
- Drying system
- Air exhaust system
- Water paths
- Water treatment
- Which process chemicals are used and how dosing works
- Blocked functions of the machine
- What are the performance limits of the machine
- Load arrangement
- Positions in the machine used for temperature measurement, positions in the machine where the temperature is reached fastest and slowest
- Process times

Furthermore, for the reprocessing of flexible endoscopes in endoscopy washer-disinfectors, the validator must be familiar with specific functions in accordance with

DIN EN ISO 15883-4, including how the endoscope fit tester works, how to monitor flow through the endoscope and how to treat water for the final rinsing process.<sup>1</sup>

Section 6 of DIN 58341 explains the scope of the validation process in accordance with DIN EN ISO 15883-1, 2 and 4 with installation (IQ), operational (OQ) and performance qualification (PQ). The scope of the test is defined in the validation plan and includes:

- Product groups and families
- Which processes are used
- Validation period
- Which process chemicals are used
- Load carriers
- Medical devices to be reprocessed with reprocessing instructions in accordance with DIN EN ISO 17664.<sup>9</sup>

The installation qualification (IQ) section describes which tests should be performed when the machine is installed and which requirements apply for process water, process chemicals, compressed air, air and water filters and the installation site itself. The operational qualification (OQ) section sets out the tests for newly installed machines and the work steps for manual cleaning and disinfection. During OQ, the monitoring functions for the door and its locking mechanisms, the rinsing arms/nozzles and the error displays must be checked. Other functional tests include the calibration/adjustment of the measurement chains, a leak test, a check of the connections of the load carriers, pressure tests and other tests specific to endoscopy washer-disinfectors in line with the manufacturer's requirements.

In terms of parameters, the water volumes for cold water, warm water and demineralised water must be checked and documented, along with the process parameters of time, temperature, pressure and dosing. The performance qualification provides evidence that the process produces clean and disinfected medical devices and that these results can be reproduced. As part of the qualification process, the suitability of the load arrangement and load carriers is documented in cleaning and disinfection performance tests. Performance qualification involves checks on the cleaning, disinfection and drying of the medical devices, as well as the final rinsing process. Based on the results of the validation, the type, scope and interval of routine checks are defined.<sup>1</sup>

Finally, Appendix A to the standard DIN 58341 defines the minimum requirements for validation reports. A table that sets out the differences between process validations for washer-disinfectors, endoscopy washer-disinfectors and the standard instructions for manual cleaning and disinfection was compiled for this purpose. The table serves as an example and an aid for the operator and validator in the assessment and documentation of the results in the validation report.<sup>1</sup>

## Conclusion

The DIN 58341 standard defines the prerequisites and required knowledge to perform validations on washer-disinfectors and endoscopy washer-disinfectors processes and manual cleaning and disinfection procedures. It can therefore serve as the basis for the expansion of the validation course for validators, e.g. as module VALI C "Requirements for the performance qualification of cleaning and disinfection processes" and for certification in accordance with §5 of the German Medical Device Ordinance.

## References

1. DIN 58341 Requirements for the validation of cleaning and disinfection processes.
2. DIN EN ISO 15883 Washer-disinfectors, parts 1/2/4/5.
3. Guidelines published by DGKH, DGSV and AKI on the validation and routine monitoring of machine-based cleaning and thermal disinfection processes for medical devices; 5th edition.
4. Guidelines published by DGKH, DGSV; DGVS, DEGEA and AKI on the validation and routine monitoring of machine-based cleaning and disinfection processes for the reprocessing of thermolabile endoscopes.
5. Guidelines published by DGKH, DGSV; AKI in partnership with VAH: Guideline on the validation of manual cleaning and chemical disinfection of medical devices.
6. MPBetreibV (German Medical Device Ordinance).
7. Recommendation from the Commission for Hospital Hygiene and Infection Prevention (KRINKO) and the German Federal Institute for Drugs and Medical Devices (BfArM): "Hygiene requirements for the reprocessing of medical devices".
8. DWM test laboratory, Dr Winfried Michels, Warburg: DIN 58341 with new requirements for cleaning performance checks – Central sterilization
9. DIN EN 17664 Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices; Beuth Verlag GmbH Berlin; 2018-04.

# The interaction between reprocessing products (specifically lumened dental instruments) and reprocessing equipment (adapters in washer-disinfectors)

Stella Nehr-Werner, Ulrike Weber

Where does the washer-disinfector begin – and where does the instrument end?

The respective intended uses of the reprocessing product and the washer-disinfector clearly indicate where one device stops and another begins. However, at intersections between the two (such as the transfer point and the adapter), it can sometimes be challenging to distinguish the wood from the trees.

Rather than focusing on the regulatory aspects of this topic from the perspective of the new Medical Device Regulation (MDR), this article looks at the practicalities of adapting medical devices to washer-disinfectors to prepare them for cleaning and disinfection.

Which intersections arise as a result of the standards for washer-disinfectors and instruments?

DIN EN ISO 17664-1 describes the information that the instrument manufacturer must provide to enable the operator to reprocess its devices properly. Part 1 encompasses not only products that can be sterilised, but also instruments that are used, directly or indirectly, on patients. Unlike Part 2, it excludes non-critical instruments. Part 1 requires the manufacturer to provide a description of at least one validated machine-based procedure, as long as the medical device can withstand this procedure. In this description, the instrument manufacturer must provide reprocessing instructions that specify the adapter/connectors to be used to achieve the correct reprocessing result.<sup>1</sup>

The DIN EN ISO 15883 series of standards is essentially a “machine standard” that sets out the minimum criteria that machines must meet for manufacturers of washer-disinfectors. Alongside the general performance requirements for washer-disinfectors and the associated accessories used for cleaning and disinfection, Part 1 also makes reference to the need to provide any necessary adapters/connectors and suitable load carriers. It also describes the procedures and measurement devices used for validation, routine monitoring and periodic testing.<sup>2</sup>

What is the reprocessing item (instrument) used for and how should it be reprocessed?

The purpose of an instrument is defined in the “intended use” of the medical device. This intended use is part of the approval; it is determined as part of the clinical assessment and is included in the operating instructions and other information provided by the manufacturer.<sup>3</sup> The intended use not only indicates the medical purpose or application of the device, but also how it is designed to be used. This defines the use of the medical device from the perspective of approval only. The risk assessment and categorisation of medical devices, as discussed in section 1.2.1 of the recommendations published by KRINKO<sup>4</sup>, are based on the origin and future location of use of the medical device, as well as on the nature and complexity of the reprocessing of the medical device. This information produces the requirements for reprocessing and, where applicable, a requirement for connectors or adapters for more complex instruments.

This is also true for the reprocessing object (washer-disinfector, adapter, components) for cleaning and disinfection. Performance is determined in the type test and general suitability is evaluated in the clinical assessment based on a number of criteria, including literature studies, validation data and cleaning attempts. Adapters and connectors may be defined as integral components of the washer-disinfector, as components of the instrument for adaptation, or as separate accessories.

Why is reprocessing dental instruments particularly challenging, given their intended use (as a medical device) and the level of contamination?


One of the challenges involved in reprocessing is the level of soiling of the instruments following treatment. In addition to organic residues such as blood, pus or saliva, anorganic materials such as cement residues, adhesives or fillers, as well as dust from abrasion and oil residues, can present a challenge during reprocessing (fig. 2). With more complex instruments and rough surfaces in particular, this soiling is often located in areas that are very difficult to access, which makes them more difficult to clean.

## Authors

Stella Nehr-Werner  
Global Clinical Affairs Manager  
Dentsply Sirona  
Fabrikstraße 31, 64625 Bensheim, Germany  
Stella.Nehr-Werner@dentsplysirona.com  
www.dentsplysirona.com

Dr Ulrike Weber  
Scientific Application & Sales Support  
Customer Segments & Solutions:  
Miele Professional business unit  
Carl-Miele-Str. 29, 33332 Gütersloh, Germany  
ulrike.weber@miele.com  
www.miele-professional.com

**Fig. 1:** Sample load arrangement.<sup>5</sup>

With lumened instruments, adapters or connectors must be used to ensure that the inner surfaces, pipes and/or hinges and joints can be thoroughly cleaned.

Anorganic soiling such as cement residues and fillers must be removed from the surface of the instruments before reprocessing and adaptation to the washer-disinfector to prevent hardening.

The process that takes place in the washer-disinfector is designed to remove organic soiling and the associated microbiological contamination; it cannot remove anorganic soiling (after all, filler is designed to remain in the mouth for many years and should not fall off when you brush your teeth).

Transmission instruments and moving parts on scissors and forceps must be treated with oil of a specific quality to ensure that they function smoothly and correctly, without wear corrosion and any loss of transmission energy. The process of reprocessing will also remove these substances, so they must be reapplied before sterilisation.

### Why is reprocessing dental instruments particularly challenging in terms of the washer-disinfector?

It is not just the washer-disinfectors themselves that have to satisfy basic standards for cleaning and disinfection, such as the DIN EN ISO 15883 series; the relevant connectors and adapters or other load carriers must also meet certain performance requirements. In an ideal scenario, there would be a finely tuned interaction between the instrument manufacturer and the device manufacturer. It is important to remember that third parties may also supply load carriers; this is often the case with special instruments.

In terms of adaptation, the focus should not only be on individual instruments – turbines, aspirators, scaler tips, multi-func-

**Fig. 2:** Dental turbine with push button removed – rotor with residual soiling.

tional syringes and similar – but also on the load and loading pattern as a whole. For reasons of simplicity, resource conservation and cost-effectiveness – and to ensure that staff are on board with the solution – the following considerations should take priority:

- Reprocessing of lumened instruments and other instruments in the same cycle in the washer-disinfector
- Option to add further adapters for lumened instruments to prepare for the possibility of practice expansion or the introduction of new special instruments (e.g. multi-functional syringe tips)
- Flexible and modular layout of baskets and inserts (e.g. fig. 1)
- Adapter and basket solutions that allow for the adaptation of various reprocessing goods (e.g. dental turbines, dental handpieces) without having to constantly change the adapter or load carrier

The procedure for reprocessing 3F cannulas and multi-functional syringe tips (MF) is a good example of how existing adapter solutions and load carriers can be combined and extended to other instruments.

In dentistry, multi-functional syringes (MF) are used for tasks such as drying and rinsing preparation sites. Generally, these syringes are permanently fixed to the treatment unit and have an air outlet and a water outlet. The tip is removed for reprocessing after each patient, and care must be taken to ensure that the tip is properly cleaned and disinfected both internally and externally using the selected reprocessing procedure (example of adaptation in fig. 4)

Adapters and connectors for lumened instruments serve as a transfer and connection point between the washer-disinfector and the reprocessing product. These components play a key role in reprocessing as they transfer (mechanical) pressure. The balance in the irrigation connectors both challenges and affects the entire process. Washer-disinfector manufacturers provide information on how a sufficiently standardised rinsing pressure can be guaranteed. Example: IFU A 105/1: [...] all screw faces must be equipped with nozzles, adapters, irrigation sleeves or blind stoppers. Damaged irrigation connectors such as nozzles, adapters or irrigation sleeves must not be used. Irrigation connectors that are not being used do not need to be replaced by blind stoppers. [...]<sup>5</sup>

The operating instructions for washer-disinfectors or adapters tell the user whether the adapters must be fitted with instru-

**Fig. 3:** Rotor of a dental turbine with rust-coloured deposits due to inadequate oiling.


ments or can have an open outlet. Modern washer-disinfectors have taken this into account during usability testing, and the operating instructions state whether and how the adapters must be closed (to ensure that pressure does not fall or rise too sharply), as well as providing instructions on loading patterns and other information (e.g. adapters may remain open as long as the thread on the basket is not visible). It is not possible to make sweeping statements such as “always close” or “can stay open”, as the technical specification (e.g. the circulation pump in the washer-disinfector) and the circulation mechanics and balance of the process in washer-disinfectors varies.

During the rinsing process, the adapter and the connected lumened instruments collect and transport liquor volume that may then no longer be available to the other elements of the rinsing mechanism (e.g. spray arm, circulation pump). However, a balanced program that is tailored to the load items is required to ensure even wash pressure for internal and external rinsing. Too little water in the washing circuit will lead to foam and a reduction in pressure.

Luer Lock adapters are the classic adapter, as the connection is standardised across the industry.

The adapter is designed not just to connect the lumened instruments, but also to ensure that the instruments can be optimally positioned (e.g. upright or at a slight angle to ensure that water can flow away easily between the rinse steps), with no undercutting and as few direct contact surfaces as possible with the instrument, as well as high levels of usability and a long service life.

If instruments require further treatment after cleaning and disinfection, a suitable adapter must also be selected for the application of the care product. For example, dental transmission instruments need to be a perfect fit for the transmission path. However, for this instrument group, there is also an option to

**Fig. 4:** Holder for reprocessible cannulas from 3F and MF syringes (Miele A 865).

use a so-called combination device, which combines cleaning, disinfection and care in one automated entire process.

## Conclusion

Ideally, the washer-disinfector and its adapters or connectors should form a single unit with the reprocessed medical devices to ensure that the end result of the process is a fully cleaned and disinfected medical device. It is always a good idea to check the operating instructions for the load and the machine to achieve the best results for the user's own typical loading patterns. The validation of the overall process with the selected combination will provide information on the efficacy and reproducibility of the process.

## References

1. DIN EN ISO 17664-1:2021-02 – draft. Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices – Part 1: Critical and semi-critical medical devices (ISO/FDIS 17664-1:2020); German and English version prEN ISO 17664-1:2021.
2. DIN EN ISO 15883-1:2021-01 – draft. Washer-disinfectors – Part 1: General requirements, terms and definitions and tests (ISO/DIS 15883-1:2020); German and English version prEN ISO 15883-1:2020.
3. REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC. Article 2 (12).
4. Commission for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute (RKI) and the German Federal Institute for Drugs and Medical Devices (BfArM): Hygiene requirements for the reprocessing of medical devices. German Federal Health Bulletin, Health research, Health protection, 2012; 55:1244-1310; pp. 1267–1269.
5. Operating instructions for injector upper basket, A 105/1. Miele material no. 11 309 120.

**Dr Hubert Holz**  
Medical specialist for  
hygiene &  
environmental medicine/  
anaesthesia

## “Three questions for ...”

### *Clostridioides difficile*

1. *What does the bacilliform bacterium Clostridioides difficile have to do with antibiotics?*

Many of the antibiotics we use – and that are often very useful in healthcare – make it easier for *Clostridioides difficile* to survive. Other important bacteria that occur naturally in the human intestinal tract are killed off, while the resistant *Clostridioides difficile* bacteria can multiply unhindered. When these bacteria dominate the gut, this can result in very serious diarrhoea.

2. *Why is Clostridioides difficile such a big problem in hospitals?*

Patients in hospital with severe bacterial infections often require antibiotic treatment. This significantly increases the aforementioned risk of *Clostridioides*-induced diarrhoea. If effective treatment to combat *Clostridioides difficile* is not administered quickly, the infection can be fatal; this has happened many times. Unfortunately, there is not always sufficient knowledge of how to effectively treat *Clostridioides difficile* infections.

3. *Where do current hygiene concepts to tackle the bacteria fall down?*

*Clostridioides difficile* is highly stable, including in uninhabited environments, primarily due to its ability to form spores. These spores are one of the very few microorganisms that we cannot combat or deplete effectively with alcohol hand disinfectants. For this bacteria, we need to wash our hands after disinfection to ensure that the spores are not transferred. When disinfecting surfaces, we also need to use special disinfectants that are effective against spores.

**Ines Korschake**  
Hygiene management

## New advisory body member: Ines Korschake

In 2004, Ines Korschake joined the Johanniter-Krankenhaus Genthin-Stendal as a state-approved Lead Hygiene Specialist, where she established a hygiene management structure. The specialist anaesthesia and intensive care nurse implemented hygiene standards based on the latest RKI specifications and heads up the team of hygiene representatives at the hospital.

Her work focuses on prevention and the promotion of patient health in the areas of hygiene, care and wound management.

Since 2012, she has been working as a lecturer at the Altmark Campus at Johanniter-Akademie Mittel-

deutschland, helping to educate the next generation of healthcare staff. She also lectures Germany's next generation of hygiene specialists at the Christliche Akademie Halle-Dörlau.

Ines places great importance on networking with other experts in our increasingly globalised world, playing an active role in the multi-resistant pathogen network of Saxony-Anhalt (HYSA) and in the working group of the University of Magdeburg. Since the outbreak of the coronavirus pandemic, she has been part of the CORONA COMPETENCE TEAM in the Johanniter organisation.

## 25th annual congress of the DGSV in Fulda, 4 to 5 October 2021

The 25th annual congress of the German Society for Sterile Supplies will this year once again focus on topics related to the reprocessing of medical products in CSSDs.

With a host of interesting presentations and discussions, participatory workshops and a parallel session for registered doctors, the congress will be a great opportunity for experts from a range of disciplines to come together and share their knowledge and experiences. Register at: [www.dgsv-kongress.de](http://www.dgsv-kongress.de)

Due to the current Covid-situation the congress organization decided to make it a virtual congress.

## 21st World Sterilization Congress in Geneva, 17 to 20 November 2021

In November 2021, the WFHSS and the Swiss Society for Sterile Supply will host the 21st World Sterilization Congress in Geneva, Switzerland. The congress will focus on sharing knowledge around the entire sterilisation process, and will include presentations on findings from the latest scientific research, new technologies and opportunities to discuss experiences and opinions.

The congress will bring together hygiene specialists and renowned experts from all over the world to discuss the latest developments and innovations from the field of sterilisation. The event is the largest international meeting focusing on hospitals, sterilisation and research.

View the programme here: <https://www.wfhss-congress.com/wfhss-2021-call-for-abstracts>

## | Legal notice

### Scientific advisory council:

H. Biering, Düsseldorf  
F. Brill, Hamburg  
J. Gebel, Bonn  
A. Hartwig, Berlin  
H. L. Holz, Mainz  
T. Miorini, Graz  
U. Junghannß, Köthen  
S. Kaufmann, Saarbrücken  
I. Korschake, Stendal  
M. Pietsch, Mainz  
B. Wilbrandt, Berlin

### Publisher:

Office, das Büro der aseptica  
Bernd Vieregge  
Frieda-Nadig-Straße 53  
33332 Gütersloh, Germany  
E-mail: [info@aseptica.com](mailto:info@aseptica.com)

### Responsible for content:

Dr Ulrike Weber  
Miele Professional business unit  
Miele & Cie. KG  
Carl-Miele-Straße 29  
33332 Gütersloh, Germany  
Tel.: +49 5241 89 1494  
E-mail: [ulrike.weber@miele.com](mailto:ulrike.weber@miele.com)

### Overall production:

COLLET Concepts Communication  
Ziethenstraße 10  
33330 Gütersloh, Germany  
Tel.: +49 5241 50 56 664  
E-mail: [info@aseptica.com](mailto:info@aseptica.com)  
Website: [www.aseptica.com](http://www.aseptica.com)  
Stefan Collet, Sandra Acikportali

### In co-operation with:

Ecolab Deutschland GmbH  
Ecolab-Allee 1 | 40789 Monheim am Rhein, Germany;  
Miele & Cie. KG  
P.O. box | 33325 Gütersloh, Germany;  
Dentsply Sirona Deutschland GmbH  
Fabrikstraße 31 I 64625 Bensheim, Germany;  
Xylem Analytics Germany Sales GmbH & Co. KG  
Ebro  
Peringerstraße 10 | 85055 Ingolstadt, Germany;  
Innovations Medical Vertriebs GmbH  
Badstraße 11 | 78532 Tuttlingen, Germany

### Editorial team:

Aaron Papadopoulos, Ecolab  
Ulrike Weber, Miele  
Stella Nehr-Werner, Dentsply Sirona  
Iven Kruse, ebro  
Michael Schändlinger, Innovations Medical

Title image: Steelco

Circulation: 6500

Publication schedule: three times a year

Printed on chlorine-free bleached paper

Only to be reprinted with the permission of the editorial team. Articles by named authors do not necessarily reflect the opinion of the editorial team. No liability is assumed for unsolicited manuscripts and photographs. The editorial team reserves the right to shorten letters from readers.

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